

**Composition**

**P08-30100 MEM Amino Acids Solution (50x)  
w/o: L-Glutamine**

	Components	mg/L
Amino acids	L-Arginine x HCl	6,320.00
	L-Cystine	1,200.00
	L-Glutamine	0.00
	L-Histidine x HCl x H <sub>2</sub> O	2,100.00
	L-Isoleucine	2,625.00
	L-Leucine	2,620.00
	L-Lysine x HCl	3,625.00
	L-Methionine	755.00
	L-Phenylalanine	1,650.00
	L-Threonine	2,380.00
	L-Tryptophan	510.00
	L-Tyrosine	1,800.00
	L-Valine	2,340.00