

Composition

**P10-028100 Trypsin 0.25%/ 1 mM EDTA 4 Na in PBS
w/o: Ca and Mg**

	Components	mg/L	mM
Inorganic Salts	Potassium chloride	200.00	2.6827632
	Potassium dihydrogen phosphate	200.00	1.4696157
	Sodium chloride	8,000.00	136.8925394
	Di-Sodium hydrogen phosphate	1,150.00	8.1008735
Other Components	EDTA-Na ₄ x 2H ₂ O	420.00	1.0091302
	Trypsin	2,500.00	0.1050420