

Zusammensetzung

| Trypsin 0,25%ig, EDTA 1 mM in HBSS, w/o Ca, Mg, w: Phenolred | | P10-029100 |
|-----------------------------------------------------------------|-----------------------------------------|------------|
| | Components | mg/L |
| Inorganic Salts | Kaliumchlorid | 400 |
| | Kaliumdihydrogenphosphat | 60 |
| | Natriumchlorid | 8000 |
| | <i>di</i> -Natriumhydrogenphosphat wfr. | 47,88 |
| | <i>D</i> (+)-Glucose wfr. | 1000 |
| | <i>Phenolrot</i> | 10 |
| Other Components | EDTA-Na ² | 416 |
| | Trypsin | 2500 |