

Composition

| | | |
|--|---|----------------------------------|
| Trypsin 0,05 %ig, EDTA 0,02 % in HBSS, w/o Ca, Mg, w: Phenolred | | P10-038100 P10-038500 mg/L |
| Components | | |
| Inorganic Salt | Kaliumchlorid | 400 |
| | Kaliumdihydrogenphosphat | 60 |
| | Natriumchlorid | 8000 |
| | <i>di</i> -Natriumhydrogenphosphat wfr. | 47,88 |
| | D(+)-Glucose wfr. | 1000 |
| | Phenolrot | 10 |
| Other Components | EDTA | 200 |
| | Trypsin | 500 |

DOMINIQUE DUTSCHER SAS