

Composition

P10-040100 **Trypsin 0.05%/ EDTA 4 Na 0.02% in HBSS**
w/o: Ca and Mg
w: Phenol red

	Components	mg/L
Inorganic Salts	Potassium chloride	400.00
	Potassium dihydrogen phosphate	60.00
	Sodium chloride	8,000.00
	Di-Sodium hydrogen phosphate anhydrous	47.88
Other Components	D(+)-Glucose dried	1,000.00
	EDTA-Na ₂	220.00
	Phenol red	10.00
	Trypsin	500.00