

SUBJECT: Addition of specifications in certificate of analysis of S0750 New Born Calf Serum

1- INFORMATION ON PRODUCT CONCERNED BY THE NOTIFICATION:

Catalogue Number: S0750

Batch Number: NA

Production date: NA

2- DESCRIPTION OF THE NOTIFICATION:

Description: In our ongoing commitment to quality and precision, we have undertaken a review to enhance the accuracy of the certificate of analysis for the product S0750. The modifications aim to provide clearer and more precise information.

Specifications for osmolality, pH and Total protein for S0750 New Born Calf Serum are currently "For information" on the certificate of analysis.

Dutscher has decided to update the specifications as follows:

- pH : 7,0 - 8,3
- Osmolality : 250 - 355 mOsm/kg
- Total Protein : 50 - 80 g/l

Documents linked: S0750 Technical Data Sheet

3- PREVENTIVE / CORRECTIVE ACTION TO REALISE:

Description: NA

Written by: E. Merlet

Position: QCEHS
Coordinator

Date: 21/12/2023

Signature:

