

SUBJECT: Update of the pH specifications of S2500 - Rabbit Serum

1- INFORMATION ON PRODUCT CONCERNED BY THE NOTIFICATION:

Catalogue Number: S2500

Batch Number: NA

Production date: Batches produced after 05/01/2024

2- DESCRIPTION OF THE NOTIFICATION:

Description:

Dutscher has analysed the pH value of all the batches of material S2500 produced since 2019. All results are inside specifications (7,50 to 8,00) but the average value is 7,78, minimum value is 7,50 and maximum value is 7,95. So the specifications are a little bit too strict.

As a consequence, Dutscher has decided to change its specifications to 7,3 - 8,3.

This change has no impact on product quality.

Documents linked: S2500 Technical Data Sheet

3- PREVENTIVE / CORRECTIVE ACTION TO REALISE:

Description: NA

Written by: E. Merlet

Position: QCEHS
Coordinator

Date: 05/01/2024

Signature:

