

Trypsin 0.05% (1x) Solution (SH30236)	
Description	
Inorganic Salts	mg/L
KCl	400.0000
KH ₂ PO ₄	60.0000
MgSO ₄ (anhydrous)	0.0000
NaCl	8000.0000
Na ₂ HPO ₄ (anhydrous)	47.6800
CaCl ₂	0.0000
Other Components	
D-Glucose	1000.0000
Phenol Red (Sodium)	11.0000
EDTA-4Na	20.0000
Trypsin (irradiated)	50.0000
NaHCO ₃	350.0000