



Formulation 29-1322-59 AA

Cell culture

HyClone™ classical media

Phosphate buffered saline solutions

Component description	Code number SH30013		Code number SH30028		Code number SH30264		Code number SH30378	
	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
Calcium chloride (anhydrous)	0	0	0	0	100	0.9011	0	0
Potassium chloride	200	2.6828	200	2.6828	200	2.6828	2000	26.8276
Potassium phosphate monobasic (anhydrous)	200	1.4696	200	1.4696	200	1.4696	2000	14.6962
Magnesium chloride-6H ₂ O	0	0	0	0	100	0.4919	0	0
Sodium chloride	8000	136.8925	8000	136.8925	8000	136.8925	80000	1368.9254
Sodium phosphate dibasic (anhydrous)	1150	8.1009	1150	8.1009	1150	8.1009	0	0
Sodium phosphate dibasic-7H ₂ O	0	0	0	0	0	0	21600	80.576

Component description	Code number SH30597		Code number SH30256		Code number SH30258	
	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
Calcium chloride (anhydrous)	1000	9.0106	0	0	0	0
Potassium chloride	2000	26.8276	0	0	0	0
Potassium phosphate monobasic (anhydrous)	2000	14.6962	144	1.0581	1440	10.5812
Magnesium chloride-6H ₂ O	1000	4.9188	0	0	0	0
Sodium chloride	80000	1368.9254	9000	154.0041	90000	1540.0411
Sodium phosphate dibasic (anhydrous)	11500	81.0087	795	5.6002	7950	56.0017
Sodium phosphate dibasic-7H ₂ O	0	0	0	0	0	0

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden
www.gelifesciences.com/hyclone

GE and GE monogram are trademarks of General Electric Company.
HyClone is a trademark of General Electric Company or one of its subsidiaries.
© 2014 General Electric Company—All rights reserved. First published Oct. 2014
All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.
GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK
GE Healthcare Europe, GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany
GE Healthcare Bio-Sciences Corp., 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA
GE Healthcare Japan Corporation, Sanken Bldg., 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan
29-1322-59 AA 10/2014