

Ivan Aston
Technical Support
HyClone Cell Culture

925 West 1800 South
Logan, Utah 84321
USA

T +1 435 792 8000
HyClone.techsupport@cytiva.com

30 September 2021

Dear Customer:

Penicillin/Streptomycin Solution, Catalog Number SV30010, contains 10000 units of Penicillin G per mL, 10000 µg Streptomycin per mL, and 8.5 g/L of NaCl (normal saline). This product is 0.1 µm sterile filtered.

Regards,



Ivan Aston
Technical Support Specialist

DOMINIQUE DUTSCHER