

**Product code : X0112**

**Product name : Earle's Balanced Salts 10x w/ Calcium w/Magnesium w/o Sodium Bicarbonate**

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	2.65000000
7487-88-9	Magnesium Sulfate Anhydrous	0.97670000
7447-40-7	Potassium Chloride	4.00000000
7647-14-5	Sodium Chloride	68.00000000
7558-80-7	Sodium Phosphate Monobasic Anhydrous	1.22000000
50-99-7	D-Glucose Anhydrous	10.00000000
34487-61-1	Phenol Red Sodium Salt	0.11000000
WATER		913.04330000