

Product code : X0113

Product name : Earle's Balanced Salts 10x w/o Calcium w/o Magnesium w/o Sodium Bicarbonate

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	4.00000000
7647-14-5	Sodium Chloride	68.00000000
7558-80-7	Sodium Phosphate Monobasic Anhydrous	1.22000000
50-99-7	D-Glucose Anhydrous	10.00000000
34487-61-1	Phenol Red Sodium Salt	0.11000000
WATER		916.67000000

DOMINIQUE DUTSCHER SAS