

Product code : X0507

Product name : HBSS 10X w/o Calcium w/o Magnesium w/o Sodium Bicarbonate w/o Phenol Red

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	4.00000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.60000000
7647-14-5	Sodium Chloride	80.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.47880000
50-99-7	D-Glucose Anhydrous	10.00000000
WATER		904.92120000