

Product code : X0509

Product name : HBSS 10X w/ Calcium w/ Magnesium w/o Sodium Bicarbonate w/ Phenol Red

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	1.85000000
7487-88-9	Magnesium Sulfate Anhydrous	0.97670000
7447-40-7	Potassium Chloride	4.00000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.60000000
7647-14-5	Sodium Chloride	80.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.47880000
50-99-7	D-Glucose Anhydrous	10.00000000
34487-61-1	Phenol Red Sodium Salt	0.11000000
WATER		901.98450000