

Product code : X0510

Product name : HBSS 10X w/o Calcium w/o Magnesium w/ Sodium Bicarbonate w/o Phenol Red

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	4.00000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.60000000
7647-14-5	Sodium Chloride	80.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.47880000
50-99-7	D-Glucose Anhydrous	10.00000000
144-55-8	Sodium Bicarbonate	3.50000000
WATER		901.42120000

DOMINIQUE DUTSCHER SAS