

**Product code : X0513**

**Product name : HBSS 10X w/o Calcium w/o Magnesium w/o Sodium Bicarbonate w/  
Phenol Red**

| CAS Number | Components                              | Quantity in g/l |
|------------|---|-----------------|
| 7447-40-7  | Potassium Chloride                      | 4.00000000      |
| 7778-77-0  | Potassium Phosphate Monobasic Anhydrous | 0.60000000      |
| 7647-14-5  | Sodium Chloride                         | 80.00000000     |
| 7558-79-4  | Sodium Phosphate Dibasic Anhydrous      | 0.47880000      |
| 50-99-7    | D-Glucose Anhydrous                     | 10.00000000     |
| 34487-61-1 | Phenol Red Sodium Salt                  | 0.11000000      |
| WATER      |   | 904.81120000    |

DOMINIQUE DUTSCHER SAS