

Product code : X0515

Product name : Dulbecco's Phosphate Buffered Saline 10x w/o Calcium w/o Magnesium

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	2.00000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	2.00000000
7647-14-5	Sodium Chloride	80.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	11.50000000
WATER		904.50000000

DOMINIQUE DUTSCHER CAS