

**Product code : X0520**

**Product name : Dulbecco's Phosphate Buffered Saline 10x w/ Calcium w/ Magnesium**

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	1.33000000
7791-18-6	Magnesium Chloride Hexahydrate	1.00000000
7447-40-7	Potassium Chloride	2.00000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	2.00000000
7647-14-5	Sodium Chloride	80.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	11.50000000
WATER		902.17000000

DOMINIQUE DUTSCHER SAS