

Product code : X0556

Product name : MEM Vitamins 100X w/o L-Glutamine

CAS Number	Components	Quantity in g/l
7647-14-5	Sodium Chloride	8.50000000
67-48-1	Choline Chloride	0.10000000
137-08-6	D-Ca Pantothenate	0.10000000
59-30-3	Folic Acid	0.10000000
87-89-8	Myo-Inositol	0.20000000
98-92-0	Nicotinamide (Nicotinic acid amide)	0.10000000
65-22-5	Pyridoxal Hydrochloride	0.10000000
83-88-5	Riboflavin	0.01000000
67-03-8	Thiamine Hydrochloride	0.10000000
WATER		990.69000000