

---

*Product code : X0557*

*Product name : MEM Non Essential Amino Acids 100Xw/o L-Glutamine*

---

CAS Number	Components	Quantity in g/l
56-40-6	Glycine	0.75000000
56-41-7	L-Alanine	0.89000000
5794-13-8	L-Asparagine Monohydrate	1.50000000
56-84-8	L-Aspartic acid	1.33000000
56-86-0	L-Glutamic Acid	1.47000000
147-85-3	L-Proline	1.15000000
56-45-1	L-Serine	1.05000000
WATER		991.86000000