

Product code : X0915

Product name : Trypsin 2.5% in PBS w/o Calcium w/o Magnesium w/o Phenol Red

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	0.20000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.20000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	1.15000000
50-99-7	D-Glucose Anhydrous	1.00000000
9002-07-7	Trypsin 1:250 (Non-ACF)	25.00000000
WATER		964.45000000

DOMINIQUE DUTSCHER CAS