

Product code : X0920

Product name : Trypsin 2.5% in HBSS w/o Calcium w/oMagnesium w/o Phenol Red

CAS Number	Components	Quantity in g/l
7447-40-7	Potassium Chloride	0.40000000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.06000000
7647-14-5	Sodium Chloride	8.00000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.04788000
50-99-7	D-Glucose Anhydrous	1.00000000
9002-07-7	Trypsin 1:250 (Non-ACF)	25.00000000
144-55-8	Sodium Bicarbonate	0.35000000
WATER		965.14212000

DOMINIQUE DUTSCHER SAS