



Oakton Conductivity Calibration Pouches, 1413µS; 20/Pk

Oakton – N° d'article FV-35653-11

Sealed for freshness and easy travel

- Tear off a corner, or the entire top—opens wide enough for all probe types
- Notched corners for easy opening from either side

Caractéristiques & Description

Conductivity Value (µS)	1413
Resistivity Value (MΩ)	0.00071
Volume (mL)	20
Volume (oz)	0.7
Volume	20 pouches
Conductivity-to-TDS Calibration Values (ppm 442)	1000
Conductivity-to-TDS Calibration Values (ppm KCl)	745
Conductivity-to-TDS Calibration Values (ppm	702

Single-use 20-mL pouches are sealed to ensure freshness—take these convenient standard solutions anywhere! Don't worry about transferring solution into another container—calibrate directly from the pouch.

A µS-to-ppm conversion table is printed on each pack for quick reference.

NaCl)	
Description	conductivity calibration pouches, 1413 μ S, box of 20 pouches
