

Catalog Number: 7286

Product Description: Corning Denatured Collagen Microcarriers, Dissolvable, 5g

Component Materials:

Beads - Polygalacturonic Acid Sodium Salt (PGA)
Surface - Collagen Type A – Porcine Gelatin, 300 Bloom
Relative Density - 1.02 ±0.01 vs H₂O

Product Dimensions:

Microcarrier Diameter - 200-300 µm Minimum weight per vial - 5g

Irradiation: This product has been irradiated at a minimum of 15.0kGy.

BSE/TSE – Product is not manufactured with materials derived from TSE relevant species as defined in EMA / 410 / 01.

Quality Control Testing – Representative production samples are collected and inspected in accordance with current applicable product specifications. Inspection records are reviewed and approved by qualified personnel for product release. Key inspections and inline tests are listed below:

Visual Inspection - Pass
Packaging Inspection - Pass

Lot Number Designation:

8 Digit Lot Number: First 3 digits – Julian date, start of manufacturing; Next 2 digits – Year of manufacture; Last 3 digits – Batch identification.

Rev No: 1