



Declaration of Conformity



Product Name : Tri-Cornered Beakers
Product Description : Biologix Tri-Cornered Beakers, 250ml, PP, Three pour spouts, Disposable, 25 Pieces/Sleeve, 4 Sleeves/Case
Catalog Number : 93-0250
Serial No. :

This is to certify that we are the supplier of above mentioned items. Beakers are made from strong, medical-grade polypropylene and are autoclavable. Tri-corner design provides three dripless pouring spouts and is available in six different sizes to cater to various needs.



Date: Sep,2025

Authorized Signature _____


Quality Control