

## PCR Strip Tubes

### I. APPLICATIONS

1. Molecular biology
2. Cell centrifugation
3. Sample storage
4. Centrifugation of precipitates

### II. FEATURES

1. Materials: Virgin polypropylene
2. Chemically resistant to alcohols and mild organic solvents
3. Transparent wall permits easy viewing of tube contents
4. Temperature range: stable from -20°C to 100°C
5. Tested to withstand centrifugation of 6,000 RCF in a fully supported rotor
6. Certified RNase & DNase Free
7. Extra transparent caps allow for RT-PCR applications

### III. TUBE DESIGN

1. Available in 0.2ml version
2. Conical bottom design fits most standard floor model and table-top centrifuges
3. Thin walls to accelerate temperature transfer
4. Clearly visible with either clear or turbid samples
5. Available in five colors: clear, pink, green, blue and yellow

### IV. CAP DESIGN

1. 0.2ml 8-Strip PCR tubes
  - a) Available in flat-top design
  - b) Available in dome-shaped design
  - c) Available in five colors: clear, pink, green, blue and yellow

### V. PACKAGING

1. 0.2ml 8-Strip PCR tubes: 125 strips/pack, 10 packs/case

### VI. DESCRIPTION

Item	Description	Packaging	Sterile
60-0088	0.2ml 8-Strip, Flat-top, Clear	125 strips/pack, 10 packs/case	No
60-0188	0.2ml 8-Strip, Flat-top, Pink	125 strips/pack, 10 packs/case	No
60-0288	0.2ml 8-Strip, Flat-top, Green	125 strips/pack, 10 packs/case	No
60-0388	0.2ml 8-Strip, Flat-top, Blue	125 strips/pack, 10 packs/case	No
60-0488	0.2ml 8-Strip, Flat-top, Yellow	125 strips/pack, 10 packs/case	No

60-0008	0.2ml 8-Strip, Dome cap, Clear	125 strips/pack, 10 packs/case	No
60-0108	0.2ml 8-Strip, Dome cap, Pink	125 strips/pack, 10 packs/case	No
60-0208	0.2ml 8-Strip, Dome cap, Green	125 strips/pack, 10 packs/case	No
60-0308	0.2ml 8-Strip, Dome cap, Blue	125 strips/pack, 10 packs/case	No
60-0408	0.2ml 8-Strip, Dome cap, Yellow	125 strips/pack, 10 packs/case	No

## VII. STERILITY

1. N/A

## VIII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


## IX. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

## X. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H	I
60-0008	Tube	20.70	3.40	O.D. 8.13	O.D. 7.11	O.D. 6.10	O.D.5.36	9.35	71.12	-
60-0088	Tube	20.70	3.40	O.D. 8.13	O.D. 7.11	O.D. 6.10	O.D.5.36	9.35	71.12	-
60-0008C	Cap	5.51	0.76	O.D. 4.42	O.D. 5.46	80.29	5.51	8.99	7.37	5.52
60-0088C	Cap	3.28	0.76	O.D. 4.42	O.D. 5.46	77.6	3.28	8.99	7.37	5.52



60-0008

