

Drosophila Vials

I. APPLICATIONS

1. Storage accessory for Drosophila related research

II. FEATURES

1. Materials: polystyrene
2. Glass-like clarity for easy sample viewing
3. Disposable

III. VIAL DESIGN

1. Available in Narrow (25 x 95mm) and Wide (28.5 x 95mm)
2. Both narrow and wide vials are available in bulk and tray packaging

IV. PACKAGING

1. Bulk, 500 pieces/case
2. Tray, 100 vials/tray, 500 pieces/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
51-0500	Narrow, 25 x 95mm vial	Bulk, 500 pieces/case	No
51-0501	Narrow, 25 x 95mm vial	Tray, 100 pieces/tray, 5 trays/case	No
51-0800	Wide, 28.5 x 95mm vial	Bulk, 500 pieces/case	No
51-0801	Wide, 28.5 x 95mm vial	Tray, 100 pieces/tray, 5 trays/case	No

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER