

## CryoKING Handheld Scanner

### I. APPLICATIONS

1. Infrared scanning items to identify barcodes and 2D barcodes accurately

### II. FEATURES

1. Support most code standards (ex. Code39, Code128, GS1 DataMatrix, GS1-128 and etc.)
2. Ergonomic handheld design
3. Could be set to multiple batch modes and turbo mode
4. Compatible with CryoKING Crown biobanking management software
5. Collect sample information rapidly and accurately
6. USB interface

### III. PACKAGING

- 1 Scanner/ Pack

### IV. DESCRIPTION

Cat #	Description	Scan Type
CKS-0101	Wired Handheld Scanner Gun	Individual Vial

### V. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

### VI. PRODUCT PHOTOS

