

Declaration of Conformity

CE, UL & CSA Standards



We declare that the product

Precision PLUS Pipettor

Is in conformity with

EU Directives:

**89/336/EEC Electromagnetic Compatibility Directive,
amended by 92/31/EEC & 93/68/EEC
73/23/EEC Low Voltage Directive**

Standards to which conformity is declared

EMI/EMC

**Emissions Testing: FCC Part 15 / EN 55011 (Class A)
Immunity Testing: - EN 61326-1
EN 61000-3-2 Harmonics / EN61000-3-3 Flicker**

Product Safety

UL/CSA/EN 61010-1

Serial Numbers (last 5 digits) in the range of: 00100 numbered sequentially by 1 to 15000

USA Contact information:

Sequence BioTech, LLC
Mr. Daniel Seguin
3 Northern Blvd., Suite A1
Amherst, NH 03031-2329
Tel: + 603 769 0963 (USA cell direct line)

The technical documentation required to demonstrate that the product meets requirements has been verified by Bureau Veritas, and is maintained by a representative of Sequence BioTech, and is available for inspection upon request (see above contact information).

Signed:

Date:5/15/2019.....

