

Technical data sheet

Bowie-Dick-Single test pack 3FSKS611102

Article No.: 3FSKS611102	Revision No.: 2
Date: 22.11.2018	Page: 1 of 1

Field of application:	Test for daily function control and quality assurance of the steam steriliser. For single use.		
Material:	Paper		
Size:	102 x 75 x 30 mm		
Sterilisation suitability:	STEAM sterilisation		
Technical Data:	Test sheet and counter sheet: sterilisation paper Inner barrier: paper Banderole: cardboard		
Print/indicator:	Print Banderole: 1-color blue Indicator test sheet: STEAM Indicator blue Sufficient penetration with steam provided its color turns from blue to complete pink.		
Standards:	Indicator system according to ISO 11140 Part 4, Indicator sheet type 2 indicator according to ISO 11140, BSI Kitemark™		
Shipping unit:	40 pcs. packed in a cardboard box		
Shelf life:	If stored under proper conditions from the date of manufacture 4 years.		
	Temperature	Relative air humidity	
	+15°C - +25°C	35 % - 50 %	
Substances:	Our product does not contain natural latex, PVC, Bisphenol A, Colophony, TSE/BSE risk materials, SVHC Substances and Phthalates as part of the formulations.		
<p>Technical data is submitted to the usual variations of operating procedures. Therefore it is not compulsory and does not exempt from validation and operating.</p>			