

Technical data sheet

Adhesive labels with steam indicator

3FSKS675130, 3FSKS675131, 3FSKS675135, 3FSKS675140, 3FSKS675145, 3FSKS675150

Article No.: 3FSKS675130, 3FSKS675131, 3FSKS675135, 3FSKS675140, 3FSKS675145, 3FSKS675150	Revision No.: 3
Date: 13.06.2019	Page: 1 of 1

Field of application: Barcode adhesive labels for transfer printers for labeling and documentation of medical devices and sterile barrier systems

Material: Duplex Vellum white
Support web: white

Size: Core: 40 mm

Article No.:	Label in mm	Reel width in mm
3FSKS675130	55 x 31	58
3FSKS675131	55 x 33	58
3FSKS675135	58 x 35	63
3FSKS675140	60 x 31	64
3FSKS675145	60 x 44	64
3FSKS675150	63 x 40	67

Sterilisation suitability: Steam sterilisation

Technical Data: Double self-adhesive
Split in front material
Cross perforation after each label
Reel winding: lable inside, foot ahead

Print/indicator: Print:1-color blue

Steam indicator pink
The indicator color changed form pink to brown

Standards: Process indicator type 1 according to ISO 11140 part 1

Shipping unit: Reel a` 1.000 pcs.
Carton-contents: 10 reels

Shelf life: If stored under proper conditions from the date of manufacture 2 years.

Temperature	Relative air humidity
+15°C - +25°C	35 % - 50 %

Substances: Our products do not contain PVC, Bisphenol A, TSE/BSE risk materials, SVHC Substances and Phthalates as part of the formulations.

Technical data is submitted to the usual variations of operating procedures. Therefore it is not compulsory and does not exempt from validation and operating.