

## Technical data sheet

### Adhesive tape “white” with STEAM Indicator 3FKLB410710, 3FKLB410712

Article No.: 3FKLB410710, 3FKLB410712	Revision No.: 1
Date: 20.03.2024	Page: 1 of 1

<b>Field of application:</b>	Packing material for sterilisation applications Suitable for crepe, non-woven and polypropylene (SMS) wrapping material.	
<b>Material:</b>	Paper	
<b>Size:</b>	Article No.	Reel with/Reel length
	3FKLB410710	19 mm / 50 m
	3FKLB410712	25 mm / 50 m
<b>Sterilisation suitability:</b>	STEAM sterilisation	
<b>Technical Data:</b>	Paper white One side coated	
<b>Print/indicator:</b>	Print: n/a Indicator: STEAM Indicator: diagonally light green indicator lines The indicator color changed from light green to dark green / black.	
<b>Standards:</b>	Process indicator type 1 according to ISO 11140 part 1	
<b>Shipping unit:</b>	Packed in a cardboard box	
	Article No.	Carton-content
	3FKLB410710	48 reels
	3FKLB410712	36 reels
<b>Shelf life:</b>	If stored under proper conditions from the date of manufacture 2 years.	
	Temperature	Relative air humidity
	+10°C - +30°C	30 % - 70 %
<b>Substances:</b>	Our product does not contain PVC, Bisphenol A, TSE/BSE risk materials, SVHC Substances and Phthalates as part of the formulations. We confirm, that the concentration at lead, mercurial, cadmium, chromate (VI) in products, packaging or packaging components does not exceed the value of 100 ppm by weight.	
<b>Technical data is submitted to the usual variations of operating procedures. Therefore it is not compulsory and does not exempt from validation and operating.</b>		