

50.701

KIT VISITOR

1 LABCOAT, 1 BOUFFANT CAP, 1 MASK, 1 PAIR OF FOOTWEAR

FEATURES & BENEFITS

- Widely used by business visitors and kitchen ones, this disposable kit ensures a new protection for each use.
- Economical and easy to use, this kit is ideal for short-term sites visit.
- Kit visitor used as environmental control object in the hygienic approach in various industries .
- Barrier against the release of hair and squams in the workplace.



TECHNICAL DESCRIPTION

Product name: Kit visitor
 Product type : Single use, non sterile
 Color : White / Blue
 Country of origin : China

Composition	Unit weight (g)	Length (cm)	Width (cm)
1 visitor smock PE with hood & 3 buttons in front	17 um	120	150
1 paper mask with 2 folds	25 g	90	20
1 PP bouffant cap with single clip	10 g	50	-
1 blue CPE shoecover with single elastic	2,7 g	41	15

INDUSTRIES

Industrial / Restaurant & catering / Community services / Hygiene

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTIC INFORMATION

OUTER CASE SPECIFICATIONS					SPECIFICATIONS CDT	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT	EAN CODE
50.701	380 x 240 x 290	5 kg	50 boxes	3 700 381 501 017	100 kits in polybag per box	3 700 381 502 465