

DATA SHEET

Gr. 0261

•DESIGNATION

- Glass microfibres filter
- Material: 100% glass fibres
- Binder: no

•APPLICATIONS

- Suspended solids in water
- Cells collection
- Waste water analysis
- Protein filtration
- Ligand binding test
- Radio-immuno test
- Drinkable water filtration

•TECHNICAL SPECIFICATIONS

| Characteristics | Values | Unities | Standards |
|------------------------|---------------|----------------|------------------|
| Weight | 52 ± 5 % | g/sqm | 150536 |
| Thickness | 0,26 | mm | |
| Filtration velocity | 100 | sec. | Herzberg |
| Pore size | 1,2 | µm | |