


TECHNICAL DATA SHEET

Part Number	038759
Description	PCR plate 96 wells



Benefits :

- ✓ PCR plate with low profile polypropylene skirt
- ✓ Uniform and thin wall
- ✓ Designed to fit ABI / Life Technologies ® 0.2ml thermal cyclers
- ✓ Can be used in many models of PCR thermal cyclers that accept regular profiles ü Colour: natural slightly frosted
- ✓ Non-sterile, free of DNases, RNases, pyrogens and metals
- ✓ Autoclavable
- ✓ Extra low adhesion: treated to reduce protein adsorption

Technical specifications :

- ✓ Volume: 0.2 ml
- ✓ Dimensions: 126.01 x 86.01 x 22.70 mm

Sales unit :

25 PCR plate 96 wells
