

Product Name: Centrifuge Tube	Catalogue No.: T408-1
This document replaces any previous version	

1. Product Description:

- **15 ml Centrifuge tube:** Disposable, non-sterile tube
- **Compatible with:**
 - T401-10 caps (sold separately)

2. Packaging:

- **Case:** 1 bag of 1000 tubes

3. Product Specifications:

- **Material:**
 - Tube: Polystyrene
- **Temperature range:** -80°C to +80°C
- **Tubes have molded graduations** from 0.25 ml to 15 ml
- **Centrifuge resistant** up to 3 000 g.

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640 and 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **BSE / TSE:** Material is BSE / TSE Free
- **LATEX:** Material is Latex Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension and volume measurements
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*