

Product Name: Pierce-It™ Closure

Catalogue No.: T407-12

This document replaces any previous version

1. Product Description:

- **Pierceable Closure for 13mm class and plastic tubes:** Non-Sterile, disposable cap
- **Available Colors:** Black, Blue, Gray, Green, Lavender, Red, White, Yellow

2. Packaging:

- **Case:** 1 bag of 1000 units

3. Product Specifications:

- **Material:**
 - Thermoplastic Elastomer
- **Temperature range:** -80°C to +121°C.
- **Autoclavable at 121°C for up to 30 minutes**

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **LATEX:** Material is Latex Free

5. Quality Assurance:

- No presence of contamination in plastic
- Visual attributes
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging and any possible contamination

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*