

Declaration of Conformity

This European Declaration of Conformity is issued under the sole responsibility of the manufacturer.

MANUFACTURER		
Name of Company	Address	SRN
Simport Scientific Inc.	2588 Bernard-Pilon, Beloeil, Quebec, J3G 4S5, Canada	CA-MF-000045478

AUTHORIZED REPRESENTATIVE			
Name of Company	Address	SRN	Telephone/email
Emergo Europe	Westervoortsedijk 60 6827 AT Arnhem The Netherlands	NL-AR-000000116	+31.70.345.8570 - phone EmergoEurope@ul.com

PRODUCT IDENTIFICATION		
Product / Trade Name	Product Code / Catalog Number	Basic UDI-DI
Histology Cassettes	See Appendix I - List of Trade Names and Product Reference Codes	667243Histology4G
Intended Purpose		Photo
Single use Tissue and biopsy containers designed for processing and embedding. Fixation, infiltration and manual sectioning of fixed or fresh tissue samples for histological analysis. They provide a writing surface for identification and for archiving as well as a support for microtome sectioning.		N/A

IVDR RISK CLASS / COMMON SPECIFICATIONS		
Device Classification		Common Specifications
Class	A	For Fixation, processing, infiltration and embedding of fixed or fresh tissue samples.
Rule	5(c)	

Simport Scientific Inc. declares that the above-mentioned products meet the provision of the following EU legislation:

- In Vitro Diagnostic Medical Device Regulation (IVDR) (EU) 2017/746

COMPANY REPRESENTATIVE: Claudine Lafond

TITLE: Vice-President

PLACE: Beloeil, Canada

SIGNATURE: 

DATE: 21/08/2025

Appendix I - List of Trade Names and Product Reference Codes

The following products are included under this Declaration of Conformity and comply with the relevant provisions of the applicable regulations:

- HISTOSETTE® Series**
- **M490** - HISTOSETTE® Tissue Cassettes (30 Degree Angle)
 - M490-2 – White
 - M490-3 – Pink
 - M490-4 – Green
 - M490-5 – Yellow
 - M490-6 – Blue
 - M490-7 – Peach
 - M490-8 – Tan
 - M490-9 – Gray
 - M490-10 – Lilac
 - M490-11 – Orange
 - M490-12 – Aqua
 - **M491** - HISTOSETTE® Biopsy Cassettes (30 Degree Angle)
 - M491-2 – White
 - M491-3 – Pink
 - M491-4 – Green
 - M491-5 – Yellow
 - M491-6 – Blue
 - M491-7 – Peach
 - M491-8 – Tan
 - M491-9 – Gray
 - M491-10 – Lilac
 - M491-11 – Orange
 - M491-12 – Aqua
 - **M498** - HISTOSETTE® Tissue Cassettes (45 Degree Angle)
 - M498-2 – White
 - M498-3 – Pink
 - M498-4 – Green
 - M498-5 – Yellow
 - M498-6 – Blue
 - M498-7 – Peach
 - M498-8 – Tan
 - M498-9 – Gray
 - M498-10 – Lilac
 - M498-11 – Orange
 - M498-12 – Aqua
 - **M499** - HISTOSETTE® Biopsy Cassettes (45 Degree Angle)
 - M499-2 – White
 - M499-3 – Pink
 - M499-4 – Green
 - M499-5 – Yellow
 - M499-6 – Blue
 - M499-7 – Peach
 - M499-8 – Tan
 - M499-9 – Gray
 - M499-10 – Lilac
 - M499-11 – Orange
 - M499-12 – Aqua
 - **M385** - HISTOSETTE® Tissue Cassettes for PRIMERA Printers in Bulk
 - M385-2 – White
 - M385-3 – Pink
 - M385-4 – Green
 - M385-5 – Yellow
 - M385-6 – Blue
 - M385-7 – Peach
 - M385-8 – Tan
 - M385-9 – Gray
 - M385-10 – Lilac
 - M385-11 – Orange
 - M385-12 – Aqua
 - **M385T** - HISTOSETTE® Tissue Cassettes for PRIMERA Printers in QUICKLOAD™ Stack (Taped)
 - M385-2T – White
 - M385-3T – Pink
 - M385-4T – Green
 - M385-5T – Yellow
 - M385-6T – Blue
 - M385-7T – Peach
 - M385-8T – Tan
 - M385-9T – Gray
 - M385-10T – Lilac
 - M385-11T – Orange
 - M385-12T – Aqua
 - **M386** - HISTOSETTE® Biopsy Cassettes for PRIMERA Printers in Bulk
 - M386-2 – White
 - M386-3 – Pink
 - M386-4 – Green
 - M386-5 – Yellow
 - M386-6 – Blue
 - M386-7 – Peach
 - M386-8 – Tan
 - M386-9 – Gray
 - M386-10 – Lilac
 - M386-11 – Orange
 - M386-12 – Aqua
 - **M386T** - HISTOSETTE® Biopsy Cassettes for PRIMERA Printers in QUICKLOAD™ Stack (Taped)
 - M386-2T – White
 - M386-3T – Pink
 - M386-4T – Green
 - M386-5T – Yellow
 - M386-6T – Blue
 - M386-7T – Peach
 - M386-8T – Tan
 - M386-9T – Gray
 - M386-10T – Lilac
 - M386-11T – Orange
 - M386-12T – Aqua
 - **M392T** - HISTOSETTE® New Smooth Printing Surface Tissue Cassettes in QUICKLOAD™ Stack
 - M392-2T – White
 - M392-3T – Pink
 - M392-4T – Green
 - M392-5T – Yellow
 - M392-6T – Blue
 - M392-7T – Peach
 - M392-8T – Tan
 - M392-9T – Gray
 - M392-10T – Lilac
 - M392-11T – Orange
 - M392-12T – Aqua
 - **M393T** - HISTOSETTE® New Smooth Printing Surface Biopsy Cassettes in QUICKLOAD™ Stack
 - M393-2T – White
 - M393-3T – Pink
 - M393-4T – Green
 - M393-5T – Yellow
 - M393-6T – Blue
 - M393-7T – Peach
 - M393-8T – Tan
 - M393-9T – Gray
 - M393-10T – Lilac
 - M393-11T – Orange
 - M393-12T – Aqua
 - **M492** - HISTOSETTE® Tissue Cassettes
 - M492-2 – White
 - M492-2A – White Laser Printer
 - M492-3 – Pink
 - M492-3A – Pink Laser Printer
 - M492-4 – Green
 - M492-4A – Green Laser Printer
 - M492-5 – Yellow
 - M492-5A – Yellow Laser Printer
 - M492-6 – Blue
 - M492-6A – Blue Laser Printer
 - M492-7 – Peach
 - M492-7A – Peach Laser Printer
 - M492-8 – Tan

- M492-8A – Tan Laser Printer
- M492-9 – Gray
- M492-9A – Gray Laser Printer
- M492-10 – Lilac
- M492-10A – Lilac Laser Printer
- M492-11 – Orange
- M492-11A – Orange Laser Printer
- M492-12 – Aqua
- M492-12A – Aqua Laser Printer
- **M492T - HISTOSETTE® Tissue Cassettes in QUICKLOAD™ Stack (Taped)**
 - M492-2T – White
 - M492-2AT – White Laser Printer
 - M492-3T – Pink
 - M492-3AT – Pink Laser Printer
 - M492-4T – Green
 - M492-4AT – Green Laser Printer
 - M492-5T – Yellow
 - M492-5AT – Yellow Laser Printer
 - M492-6T – Blue
 - M492-6AT – Blue Laser Printer
 - M492-7T – Peach
 - M492-7AT – Peach Laser Printer
 - M492-8T – Tan
 - M492-8AT – Tan Laser Printer
 - M492-9T – Gray
 - M492-9AT – Gray Laser Printer
 - M492-10T – Lilac
 - M492-10AT – Lilac Laser Printer
 - M492-11T – Orange
 - M492-11AT – Orange Laser Printer
 - M492-12T – Aqua
 - M492-12AT – Aqua Laser Printer
- **M492-PC - HISTOSETTE® Priority Tissue Cassettes**
 - M492-16PC – Fluo Pink
 - M492-17PC – Fluo Green
 - M492-18PC – Fluo Yellow
- **M493 - HISTOSETTE® Biopsy Cassettes**
 - M493-2 – White
 - M493-2A – White Laser Printer
 - M493-3 – Pink
 - M493-3A – Pink Laser Printer
 - M493-4 – Green
 - M493-4A – Green Laser Printer
 - M493-5 – Yellow
 - M493-5A – Yellow Laser Printer
 - M493-6 – Blue
 - M493-6A – Blue Laser Printer
 - M493-7 – Peach
 - M493-7A – Peach Laser Printer
 - M493-8 – Tan
 - M493-8A – Tan Laser Printer
 - M493-9 – Gray
 - M493-9A – Gray Laser Printer
 - M493-10 – Lilac
 - M493-10A – Lilac Laser Printer
 - M493-11 – Orange
 - M493-11A – Orange Laser Printer
 - M493-12 – Aqua
 - M493-12A – Aqua Laser Printer
- **M493T - HISTOSETTE® Biopsy Cassettes in QUICKLOAD™ Stack (Taped)**
 - M493-2T – White
 - M493-2AT – White Laser Printer
 - M493-3T – Pink
 - M493-3AT – Pink Laser Printer
 - M493-4T – Green
 - M493-4AT – Green Laser Printer
 - M493-5T – Yellow
 - M493-5AT – Yellow Laser Printer
 - M493-6T – Blue
 - M493-6AT – Blue Laser Printer
 - M493-7T – Peach
 - M493-7AT – Peach Laser Printer
 - M493-8T – Tan
- **M493-PC - HISTOSETTE® Priority Biopsy Cassettes**
 - M493-16PC – Fluo Pink
 - M493-17PC – Fluo Green
 - M493-18PC – Fluo Yellow
- **M485 - HISTOSETTE® Tissue Cassettes**
 - M485-2 – White
 - M485-3 – Pink
 - M485-4 – Green
 - M485-5 – Yellow
 - M485-6 – Blue
 - M485-7 – Peach
 - M485-8 – Tan
 - M485-9 – Gray
 - M485-10 – Lilac
 - M485-11 – Orange
 - M485-12 – Aqua
- **M485SSL - HISTOSETTE® Tissue Cassettes in QUICKLOAD™ Sleeves**
 - M485-2SL – White
 - M485-3SL – Pink
 - M485-4SL – Green
 - M485-5SL – Yellow
 - M485-6SL – Blue
 - M485-7SL – Peach
 - M485-8SL – Tan
 - M485-9SL – Gray
 - M485-10SL – Lilac
 - M485-11SL – Orange
 - M485-12SL – Aqua
- **M486 - HISTOSETTE® Biopsy Cassettes**
 - M486-2 – White
 - M486-3 – Pink
 - M486-4 – Green
 - M486-5 – Yellow
 - M486-6 – Blue
 - M486-7 – Peach
 - M486-8 – Tan

- M486-9 – Gray
- M486-10 – Lilac
- M486-11 – Orange
- M486-12 – Aqua
- **M486SL - HISTOSETTE®**
Biopsy Cassettes in QUICKLOAD™ Sleeves
 - M486-2SL – White
 - M486-3SL – Pink
 - M486-4SL – Green
 - M486-5SL – Yellow
 - M486-6SL – Blue
 - M486-7SL – Peach
 - M486-8SL – Tan
 - M486-9SL – Gray
 - M486-10SL – Lilac
 - M486-11SL – Orange
 - M486-12SL – Aqua
- **M482-LI - HISTOSETTE®** - Lids Only for Tissue Cassettes in E-Z LOAD™
 - M482-2LI – White
 - M482-3LI – Pink
 - M482-4LI – Green
 - M482-5LI – Yellow
 - M482-6LI – Blue
 - M482-7LI – Peach
 - M482-8LI – Tan
 - M482-9LI – Gray
 - M482-10LI – Lilac
 - M482-11LI – Orange
 - M482-12LI – Aqua
- **M482-BA - HISTOSETTE®** - Base Only for Tissue Cassettes in E-Z LOAD™ Stack
 - M482-2BA – White
 - M482-3BA – Pink
 - M482-4BA – Green
 - M482-5BA – Yellow
 - M482-6BA – Blue
 - M482-7BA – Peach
 - M482-8BA – Tan
 - M482-9BA – Gray
 - M482-10BA – Lilac
 - M482-11BA – Orange
 - M482-12BA – Aqua
- **M483-BA - HISTOSETTE®** - Base Only for Biopsy Cassettes in E-Z LOAD™ Stack
 - M483-2BA – White
 - M483-3BA – Pink
 - M483-4BA – Green
 - M483-5BA – Yellow
 - M483-6BA – Blue
 - M483-7BA – Peach
 - M483-8BA – Tan
 - M483-9BA – Gray
 - M483-10BA – Lilac
- M483-11BA – Orange
- M483-12BA – Aqua
- **M483-LI - HISTOSETTE®** - Lids Only for Biopsy Cassettes in E-Z LOAD™ Stack
 - M483-2LI – White
 - M483-3LI – Pink
 - M483-4LI – Green
 - M483-5LI – Yellow
 - M483-6LI – Blue
 - M483-7LI – Peach
 - M483-8LI – Tan
 - M483-9LI – Gray
 - M483-10LI – Lilac
 - M483-11LI – Orange
 - M483-12LI – Aqua
- MICROMESH™ Series**
- **M407 - MICROMESH™** Biopsy Cassettes for PRIMERA Printers in Bulk
 - M407-2 – White
 - M407-3 – Pink
 - M407-4 – Green
 - M407-5 – Yellow
 - M407-6 – Blue
 - M407-7 – Peach
 - M407-8 – Tan
 - M407-9 – Gray
 - M407-10 – Lilac
 - M407-11 – Orange
 - M407-12 – Aqua
- **M407T - MICROMESH™** Biopsy Cassettes for PRIMERA Printers in QUICKLOAD™ Stack (Taped)
 - M407-2T – White
 - M407-3T – Pink
 - M407-4T – Green
 - M407-5T – Yellow
 - M407-6T – Blue
 - M407-7T – Peach
 - M407-8T – Tan
 - M407-9T – Gray
 - M407-10T – Lilac
 - M407-11T – Orange
 - M407-12T – Aqua
- **M507 - MICROMESH™** Biopsy Cassettes
 - M507-2 – White
 - M507-2A – White Laser Printer
 - M507-3 – Pink
 - M507-3A – Pink Laser Printer
 - M507-4 – Green
 - M507-4A – Green Laser Printer
 - M507-5 – Yellow
 - M507-5A – Yellow Laser Printer
 - M507-6 – Blue
 - M507-6A – Blue Laser Printer
 - M507-7 – Peach
- M507-5 – Yellow
- M507-5A – Yellow Laser Printer
- M507-6 – Blue
- M507-6A – Blue Laser Printer
- M507-7 – Peach
- M507-7A – Peach Laser Printer
- M507-8 – Tan
- M507-8A – Tan Laser Printer
- M507-9 – Gray
- M507-9A – Gray Laser Printer
- M507-10 – Lilac
- M507-10A – Lilac Laser Printer
- M507-11 – Orange
- M507-11A – Orange Laser Printer
- M507-12 – Aqua
- M507-12A – Aqua Laser Printer
- **M507SL - MICROMESH™** Biopsy Cassettes in QUICKLOAD™ Sleeves
 - M507-2SL – White
 - M507-3SL – Pink
 - M507-4SL – Green
 - M507-5SL – Yellow
 - M507-6SL – Blue
 - M507-7SL – Peach
 - M507-8SL – Tan
 - M507-9SL – Gray
 - M507-10SL – Lilac
 - M507-11SL – Orange
 - M507-12SL – Aqua
- **M507T - MICROMESH™** Biopsy Cassettes in QUICKLOAD™ Stack (Taped)
 - M507-2T – White
 - M507-2AT – White Laser Printer
 - M507-3T – Pink
 - M507-3AT – Pink Laser Printer
 - M507-4T – Green
 - M507-4AT – Green Laser Printer
 - M507-5T – Yellow
 - M507-5AT – Yellow Laser Printer
 - M507-6T – Blue
 - M507-6AT – Blue Laser Printer
 - M507-7T – Peach

- M507-7AT – Peach Laser Printer
 - M507-8T – Tan
 - M507-8AT – Tan Laser Printer
 - M507-9T – Gray
 - M507-9AT – Gray Laser Printer
 - M507-10T – Lilac
 - M507-10AT – Lilac Laser Printer
 - M507-11T – Orange
 - M507-11AT – Orange Laser Printer
 - M507-12T – Aqua
 - M507-12AT – Aqua Laser Printer
 - **M508 - MICROMESH™ Biopsy Cassettes with 4 Compartments**
 - M508-2 – White
 - M508-3 – Pink
 - M508-4 – Green
 - M508-5 – Yellow
 - M508-6 – Blue
 - M508-7 – Peach
 - M508-8 – Tan
 - M508-9 – Gray
 - M508-10 – Lilac
 - M508-11 – Orange
 - M508-12 – Aqua
- SLIMSETTE™ Series**
- **M409 - SLIMSETTE™ Tissue Cassettes for PRIMERA Printers in Bulk**
 - M409-2 – White
 - M409-3 – Pink
 - M409-4 – Green
 - M409-5 – Yellow
 - M409-6 – Blue
 - M409-7 – Peach
 - M409-8 – Tan
 - M409-9 – Gray
 - M409-10 – Lilac
 - M409-11 – Orange
 - M409-12 – Aqua
 - **M409T - SLIMSETTE™ Tissue Cassettes for PRIMERA Printers in QUICKLOAD™ Stack (Taped)**
 - M409-2T – White
 - M409-3T – Pink
 - M409-4T – Green
 - M409-5T – Yellow
 - M409-6T – Blue
 - M409-7T – Peach
 - M409-8T – Tan
 - M409-9T – Gray
 - M409-10T – Lilac
 - M409-11T – Orange
 - M409-12T – Aqua
- M409-10T – Lilac
 - M409-11T – Orange
 - M409-12T – Aqua
 - **M410 - SLIMSETTE™ Biopsy Cassettes for PRIMERA Printers in Bulk**
 - M410-2 – White
 - M410-3 – Pink
 - M410-4 – Green
 - M410-5 – Yellow
 - M410-6 – Blue
 - M410-7 – Peach
 - M410-8 – Tan
 - M410-9 – Gray
 - M410-10 – Lilac
 - M410-11 – Orange
 - M410-12 – Aqua
 - **M410T - SLIMSETTE™ Biopsy Cassettes for PRIMERA Printers in QUICKLOAD™ Stack (Taped)**
 - M410-2T – White
 - M410-3T – Pink
 - M410-4T – Green
 - M410-5T – Yellow
 - M410-6T – Blue
 - M410-7T – Peach
 - M410-8T – Tan
 - M410-9T – Gray
 - M410-10T – Lilac
 - M410-11T – Orange
 - M410-12T – Aqua
 - **M509 - SLIMSETTE™ Tissue Cassettes:**
 - M509-2 – White
 - M509-2A – White Laser Printer
 - M509-3 – Pink
 - M509-3A – Pink Laser Printer
 - M509-4 – Green
 - M509-4A – Green Laser Printer
 - M509-5 – Yellow
 - M509-5A – Yellow Laser Printer
 - M509-6 – Blue
 - M509-6A – Blue Laser Printer
 - M509-7 – Peach
 - M509-7A – Peach Laser Printer
 - M509-8 – Tan
 - M509-8A – Tan Laser Printer
 - M509-9 – Gray
 - M509-9A – Gray Laser Printer
 - M509-10 – Lilac
 - M509-10A – Lilac Laser Printer
 - M509-11 – Orange
 - M509-11A – Orange Laser Printer
 - M509-12 – Aqua
 - M509-12A – Aqua Laser Printer
 - **M509SL - SLIMSETTE™ Tissue Cassettes in QUICKLOAD™ Sleeves**
 - M509-2SL – White
 - M509-3SL – Pink
 - M509-4SL – Green
 - M509-5SL – Yellow
 - M509-6SL – Blue
 - M509-7SL – Peach
 - M509-8SL – Tan
 - M509-9SL – Gray
 - M509-10SL – Lilac
 - M509-11SL – Orange
 - M509-12SL – Aqua
 - **M509T - SLIMSETTE™ Tissue Cassettes in QUICKLOAD™ Stack (Taped)**
 - M509-2T – White
 - M509-2AT – White Laser Printer
 - M509-3T – Pink
 - M509-3AT – Pink Laser Printer
 - M509-4T – Green
 - M509-4AT – Green Laser Printer
 - M509-5T – Yellow
 - M509-5AT – Yellow Laser Printer
 - M509-6T – Blue
 - M509-6AT – Blue Laser Printer
 - M509-7T – Peach
 - M509-7AT – Peach Laser Printer
 - M509-8T – Tan
 - M509-8AT – Tan Laser Printer
 - M509-9T – Gray
 - M509-9AT – Gray Laser Printer
 - M509-10T – Lilac
 - M509-10AT – Lilac Laser Printer
 - M509-11T – Orange
 - M509-11AT – Orange Laser Printer
 - M509-12T – Aqua

- M509-12AT – Aqua Laser Printer
- **M510 - SLIMSETTE™ Biopsy Cassettes**
 - M510-2 – White Laser Printer
 - M510-2A – White Laser Printer
 - M510-3 – Pink Laser Printer
 - M510-3A – Pink Laser Printer
 - M510-4 – Green Laser Printer
 - M510-4A – Green Laser Printer
 - M510-5 – Yellow Laser Printer
 - M510-5A – Yellow Laser Printer
 - M510-6 – Blue Laser Printer
 - M510-6A – Blue Laser Printer
 - M510-7 – Peach Laser Printer
 - M510-7A – Peach Laser Printer
 - M510-8 – Tan Laser Printer
 - M510-8A – Tan Laser Printer
 - M510-9 – Gray Laser Printer
 - M510-9A – Gray Laser Printer
 - M510-10 – Lilac Laser Printer
 - M510-10A – Lilac Laser Printer
 - M510-11 – Orange Laser Printer
 - M510-11A – Orange Laser Printer
 - M510-12 – Aqua Laser Printer
 - M510-12A – Aqua Laser Printer
- **M510SL - SLIMSETTE™ Biopsy Cassettes in QUICKLOAD™ Sleeves**
 - M510-2SL – White
 - M510-3SL – Pink
 - M510-4SL – Green
 - M510-5SL – Yellow
 - M510-6SL – Blue
 - M510-7SL – Peach
 - M510-8SL – Tan
 - M510-9SL – Gray
 - M510-10SL – Lilac
 - M510-11SL – Orange
 - M510-12SL – Aqua
- **M510T - SLIMSETTE™ Biopsy Cassettes in QUICKLOAD™ Stack (Taped)**
 - M510-2T – White Laser Printer
 - M510-2AT – White Laser Printer
 - M510-3T – Pink Laser Printer
- M510-3AT – Pink Laser Printer
- M510-4T – Green Laser Printer
- M510-4AT – Green Laser Printer
- M510-5T – Yellow Laser Printer
- M510-5AT – Yellow Laser Printer
- M510-6T – Blue Laser Printer
- M510-6AT – Blue Laser Printer
- M510-7T – Peach Laser Printer
- M510-7AT – Peach Laser Printer
- M510-8T – Tan Laser Printer
- M510-8AT – Tan Laser Printer
- M510-9T – Gray Laser Printer
- M510-9AT – Gray Laser Printer
- M510-10T – Lilac Laser Printer
- M510-10AT – Lilac Laser Printer
- M510-11T – Orange Laser Printer
- M510-11AT – Orange Laser Printer
- M510-12T – Aqua Laser Printer
- M510-12AT – Aqua Laser Printer
- **M511 - SLIMSETTE™ Tissue Cassettes with 4 Compartments**
 - M511-2 – White
 - M511-3 – Pink
 - M511-4 – Green
 - M511-5 – Yellow
 - M511-6 – Blue
 - M511-7 – Peach
 - M511-8 – Tan
 - M511-9 – Gray
 - M511-10 – Lilac
 - M511-11 – Orange
 - M511-12 – Aqua
- **UNISETTE™ Series**
 - **M405 - UNISETTE™ Tissue Cassettes for PRIMERA Printers in Bulk**
 - M405-2 – White
 - M405-3 – Pink
 - M405-4 – Green
 - M405-5 – Yellow
 - M405-6 – Blue
 - M405-7 – Peach
 - M405-8 – Tan
 - M405-9 – Gray
 - M405-10 – Lilac
 - M405-11 – Orange
 - M405-12 – Aqua
 - **M405T - UNISETTE™ Tissue Cassettes for PRIMERA Printers in QUICKLOAD™ Stack (Taped)**
 - M405-2T – White
 - M405-3T – Pink
 - M405-4T – Green
 - M405-5T – Yellow
 - M405-6T – Blue
 - M405-7T – Peach
 - M405-8T – Tan
 - M405-9T – Gray
 - M405-10T – Lilac
 - M405-11T – Orange
 - M405-12T – Aqua
 - **M406 - UNISETTE™ Biopsy Cassettes for PRIMERA Printers in Bulk**
 - M406-2 – White
 - M406-3 – Pink
 - M406-4 – Green
 - M406-5 – Yellow
 - M406-6 – Blue
 - M406-7 – Peach
 - M406-8 – Tan
 - M406-9 – Gray
 - M406-10 – Lilac
 - M406-11 – Orange
 - M406-12 – Aqua
 - **M406T - UNISETTE™ Biopsy Cassettes for PRIMERA Printers in QUICKLOAD™ Stack (Taped)**
 - M406-2T – White
 - M406-3T – Pink
 - M406-4T – Green
 - M406-5T – Yellow
 - M406-6T – Blue
 - M406-7T – Peach
 - M406-8T – Tan
 - M406-9T – Gray
 - M406-10T – Lilac
 - M406-11T – Orange
 - M406-12T – Aqua
 - **M505 - UNISETTE™ Tissue Cassettes**
 - M505-2 – White
 - M505-3 – Pink
 - M505-4 – Green
 - M505-5 – Yellow
 - M505-6 – Blue
 - M505-7 – Peach
 - M505-8 – Tan
 - M505-9 – Gray
 - M505-10 – Lilac
 - M505-11 – Orange
 - M505-12 – Aqua

- **M505T** - UNISLETTE™ Tissue Cassettes in QUICKLOAD™ Stack (Taped)
 - M505-2T – White
 - M505-3T – Pink
 - M505-4T – Green
 - M505-5T – Yellow
 - M505-6T – Blue
 - M505-7T – Peach
 - M505-8T – Tan
 - M505-9T – Gray
 - M505-10T – Lilac
 - M505-11T – Orange
 - M505-12T – Aqua
- **M506** - UNISLETTE™ Biopsy Cassettes
 - M506-2 – White
 - M506-3 – Pink
 - M506-4 – Green
 - M506-5 – Yellow
 - M506-6 – Blue
 - M506-7 – Peach
 - M506-8 – Tan
 - M506-9 – Gray
 - M506-10 – Lilac
 - M506-11 – Orange
 - M506-12 – Aqua
- **M506T** - UNISLETTE™ Biopsy Cassettes in QUICKLOAD™ Stack (Taped)
 - M506-2T – White
 - M506-3T – Pink
 - M506-4T – Green
 - M506-5T – Yellow
 - M506-6T – Blue
 - M506-7T – Peach
 - M506-8T – Tan
 - M506-9T – Gray
 - M506-10T – Lilac
 - M506-11T – Orange
 - M506-12T – Aqua
- **M525T** - UNISLETTE™ Tissue Cassettes in QUICKLOAD™ Stack (Taped)
 - M525-2T – White
 - M525-3T – Pink
 - M525-4T – Green
 - M525-5T – Yellow
 - M525-6T – Blue
 - M525-7T – Peach
 - M525-8T – Tan
 - M525-9T – Gray
 - M525-10T – Lilac
 - M525-11T – Orange
 - M525-12T – Aqua
- **M526T** - UNISLETTE™ Biopsy Cassettes in QUICKLOAD™ Stack (Taped)
 - M526-2T – White
 - M526-3T – Pink
 - M526-4T – Green
 - M526-5T – Yellow
 - M526-6T – Blue
 - M526-7T – Peach
 - M526-8T – Tan
 - M526-9T – Gray
 - M526-10T – Lilac
 - M526-11T – Orange
 - M526-12T – Aqua
- **SUPERFLO™ Series**
 - **M430T** - SUPERFLO™ Tissue Cassette in QUICKLOAD™ Stack (Taped)
 - M430-2T – White
 - M430-3T – Pink
 - M430-4T – Green
 - M430-5T – Yellow
 - M430-6T – Blue
 - M430-7T – Peach
 - M430-8T – Tan
 - M430-9T – Gray
 - M430-10T – Lilac
 - M430-11T – Orange
 - M430-12T – Aqua
 - **M431T** - SUPERFLO™ Biopsy Cassette in QUICKLOAD™ Stack (Taped)
 - M431-2T – White
 - M431-3T – Pink
 - M431-4T – Green
 - M431-5T – Yellow
 - M431-6T – Blue
 - M431-7T – Peach
 - M431-8T – Tan
 - M431-9T – Gray
 - M431-10T – Lilac
 - M431-11T – Orange
 - M431-12T – Aqua
 - **M432T** - SUPERFLO™ Micro-Biopsy Cassette in QUICKLOAD™ Stack (Taped)
 - M432-2T – White
 - M432-3T – Pink
 - M432-4T – Green
 - M432-5T – Yellow
 - M432-6T – Blue
 - M432-7T – Peach
 - M432-8T – Tan
 - M432-9T – Gray
 - M432-10T – Lilac
 - M432-11T – Orange
 - M432-12T – Aqua
- **M433T** - SUPERFLO™ Tissue Cassette in QUICKLOAD™ Stack (Taped)
 - M433-2T – White
 - M433-3T – Pink
 - M433-4T – Green
 - M433-5T – Yellow
 - M433-6T – Blue
 - M433-7T – Peach
 - M433-8T – Tan
 - M433-9T – Gray
 - M433-10T – Lilac
 - M433-11T – Orange
 - M433-12T – Aqua
- **M434T** - SUPERFLO™ Biopsy Cassette in QUICKLOAD™ Stack (Taped)
 - M434-2T – White
 - M434-3T – Pink
 - M434-4T – Green
 - M434-5T – Yellow
 - M434-6T – Blue
 - M434-7T – Peach
 - M434-8T – Tan
 - M434-9T – Gray
 - M434-10T – Lilac
 - M434-11T – Orange
 - M434-12T – Aqua
- **M435T** - SUPERFLO™ Micro-Biopsy Cassette in QUICKLOAD™ Stack (Taped)
 - M435-2T – White
 - M435-3T – Pink
 - M435-4T – Green
 - M435-5T – Yellow
 - M435-6T – Blue
 - M435-7T – Peach
 - M435-8T – Tan
 - M435-9T – Gray
 - M435-10T – Lilac
 - M435-11T – Orange
 - M435-12T – Aqua
- **M530** - SUPERFLO™ Tissue Cassette
 - M530-2 – White Laser Printer
 - M530-2A – White Laser Printer
 - M530-3 – Pink Laser Printer
 - M530-3A – Pink Laser Printer
 - M530-4 – Green Laser Printer
 - M530-4A – Green Laser Printer
 - M530-5 – Yellow Laser Printer
 - M530-5A – Yellow Laser Printer
 - M530-6 – Blue Laser Printer

- M530-6A – Blue Laser Printer
- M530-7 – Peach Laser Printer
- M530-8 – Tan Laser Printer
- M530-9 – Gray Laser Printer
- M530-10 – Lilac Laser Printer
- M530-11 – Orange Laser Printer
- M530-12 – Aqua Laser Printer
- **M530AT** - SUPERFLO™ Tissue Cassette with Additive in QUICKLOAD™ Stack (Taped)
 - M530-2AT – White Laser Printer
 - M530-3AT – Pink Laser Printer
 - M530-4AT – Green Laser Printer
 - M530-5AT – Yellow Laser Printer
 - M530-6AT – Blue Laser Printer
 - M530-7AT – Peach Laser Printer
 - M530-8AT – Tan Laser Printer
 - M530-9AT – Gray Laser Printer
 - M530-10AT – Lilac Laser Printer
 - M530-11AT – Orange Laser Printer
 - M530-12AT – Aqua Laser Printer
- **M531** - SUPERFLO™ Biopsy Cassette
 - M531-2 – White Laser Printer
 - M531-3 – Pink Laser Printer
 - M531-4 – Green Laser Printer
 - M531-5 – Yellow Laser Printer
- M531-5A – Yellow Laser Printer
- M531-6 – Blue Laser Printer
- M531-7 – Peach Laser Printer
- M531-8 – Tan Laser Printer
- M531-9 – Gray Laser Printer
- M531-10 – Lilac Laser Printer
- M531-11 – Orange Laser Printer
- M531-12 – Aqua Laser Printer
- **M531AT** - SUPERFLO™ Biopsy Cassette with Additive in QUICKLOAD™ Stack (Taped)
 - M531-2AT – White Laser Printer
 - M531-3AT – Pink Laser Printer
 - M531-4AT – Green Laser Printer
 - M531-5AT – Yellow Laser Printer
 - M531-6AT – Blue Laser Printer
 - M531-7AT – Peach Laser Printer
 - M531-8AT – Tan Laser Printer
 - M531-9AT – Gray Laser Printer
 - M531-10AT – Lilac Laser Printer
 - M531-11AT – Orange Laser Printer
 - M531-12AT – Aqua Laser Printer
- **M532** - SUPERFLO™ Micro-Biopsy Cassette
 - M532-2 – White Laser Printer
 - M532-3 – Pink Laser Printer
 - M532-4 – Green Laser Printer
- M532-4A – Green Laser Printer
- M532-5 – Yellow Laser Printer
- M532-5A – Yellow Laser Printer
- M532-6 – Blue Laser Printer
- M532-6A – Blue Laser Printer
- M532-7 – Peach Laser Printer
- M532-7A – Peach Laser Printer
- M532-8 – Tan Laser Printer
- M532-8A – Tan Laser Printer
- M532-9 – Gray Laser Printer
- M532-9A – Gray Laser Printer
- M532-10 – Lilac Laser Printer
- M532-10A – Lilac Laser Printer
- M532-11 – Orange Laser Printer
- M532-11A – Orange Laser Printer
- M532-12 – Aqua Laser Printer
- M532-12A – Aqua Laser Printer
- **M532AT** - SUPERFLO™ Micro-Biopsy Cassette with Additive in QUICKLOAD™ Stack (Taped)
 - M532-2AT – White Laser Printer
 - M532-3AT – Pink Laser Printer
 - M532-4AT – Green Laser Printer
 - M532-5AT – Yellow Laser Printer
 - M532-6AT – Blue Laser Printer
 - M532-7AT – Peach Laser Printer
 - M532-8AT – Tan Laser Printer
 - M532-9AT – Gray Laser Printer
 - M532-10AT – Lilac Laser Printer
 - M532-11AT – Orange Laser Printer
 - M532-12AT – Aqua Laser Printer
- **M533** - SUPERFLO™ Tissue Cassette
 - M533-2 – White Laser Printer
 - M533-3 – Pink Laser Printer

- M533-3A – Pink Laser Printer
- M533-4 – Green
- M533-4A – Green Laser Printer
- M533-5 – Yellow
- M533-5A – Yellow Laser Printer
- M533-6 – Blue
- M533-6A – Blue Laser Printer
- M533-7 – Peach
- M533-7A – Peach Laser Printer
- M533-8 – Tan
- M533-8A – Tan Laser Printer
- M533-9 – Gray
- M533-9A – Gray Laser Printer
- M533-10 – Lilac
- M533-10A – Lilac Laser Printer
- M533-11 – Orange
- M533-11A – Orange Laser Printer
- M533-12 – Aqua
- M533-12A – Aqua Laser Printer
- **M533SL - SUPERFLO™ Tissue Cassette in QUICKLOAD™ Sleeves**
 - M533-2SL – White
 - M533-3SL – Pink
 - M533-4SL – Green
 - M533-5SL – Yellow
 - M533-6SL – Blue
 - M533-7SL – Peach
 - M533-8SL – Tan
 - M533-9SL – Gray
 - M533-10SL – Lilac
 - M533-11SL – Orange
 - M533-12SL – Aqua
- **M533AT - SUPERFLO™ Tissue Cassette with Additive in QUICKLOAD™ Stack (Taped)**
 - M533-2AT – White Laser Printer
 - M533-3AT – Pink Laser Printer
 - M533-4AT – Green Laser Printer
 - M533-5AT – Yellow Laser Printer
 - M533-6AT – Blue Laser Printer
 - M533-7AT – Peach Laser Printer
- M533-8AT – Tan Laser Printer
- M533-9AT – Gray Laser Printer
- M533-10AT – Lilac Laser Printer
- M533-11AT – Orange Laser Printer
- M533-12AT – Aqua Laser Printer
- **M534 - SUPERFLO™ Biopsy Cassette**
 - M534-2 – White
 - M534-2A – White Laser Printer
 - M534-3 – Pink
 - M534-3A – Pink Laser Printer
 - M534-4 – Green
 - M534-4A – Green Laser Printer
 - M534-5 – Yellow
 - M534-5A – Yellow Laser Printer
 - M534-6 – Blue
 - M534-6A – Blue Laser Printer
 - M534-7 – Peach
 - M534-7A – Peach Laser Printer
 - M534-8 – Tan
 - M534-8A – Tan Laser Printer
 - M534-9 – Gray
 - M534-9A – Gray Laser Printer
 - M534-10 – Lilac
 - M534-10A – Lilac Laser Printer
 - M534-11 – Orange
 - M534-11A – Orange Laser Printer
 - M534-12 – Aqua
 - M534-12A – Aqua Laser Printer
- **M534SL - SUPERFLO™ Biopsy Cassette in QUICKLOAD™ Sleeves**
 - M534-2SL – White
 - M534-3SL – Pink
 - M534-4SL – Green
 - M534-5SL – Yellow
 - M534-6SL – Blue
 - M534-7SL – Peach
 - M534-8SL – Tan
 - M534-9SL – Gray
 - M534-10SL – Lilac
 - M534-11SL – Orange
- M534-12SL – Aqua
- **M534AT - SUPERFLO™ Biopsy Cassette with Additive in QUICKLOAD™ Stack (Taped)**
 - M534-2AT – White Laser Printer
 - M534-3AT – Pink Laser Printer
 - M534-4AT – Green Laser Printer
 - M534-5AT – Yellow Laser Printer
 - M534-6AT – Blue Laser Printer
 - M534-7AT – Peach Laser Printer
 - M534-8AT – Tan Laser Printer
 - M534-9AT – Gray Laser Printer
 - M534-10AT – Lilac Laser Printer
 - M534-11AT – Orange Laser Printer
 - M534-12AT – Aqua Laser Printer
- **M535 - SUPERFLO™ Micro-Biopsy Cassette**
 - M535-2 – White
 - M535-2A – White Laser Printer
 - M535-3 – Pink
 - M535-3A – Pink Laser Printer
 - M535-4 – Green
 - M535-4A – Green Laser Printer
 - M535-5 – Yellow
 - M535-5A – Yellow Laser Printer
 - M535-6 – Blue
 - M535-6A – Blue Laser Printer
 - M535-7 – Peach
 - M535-7A – Peach Laser Printer
 - M535-8 – Tan
 - M535-8A – Tan Laser Printer
 - M535-9 – Gray
 - M535-9A – Gray Laser Printer
 - M535-10 – Lilac
 - M535-10A – Lilac Laser Printer
 - M535-11 – Orange
 - M535-11A – Orange Laser Printer

- M535-12 – Aqua
 - M535-12A – Aqua Laser Printer
 - **M535SL** - SUPERFLO™ Micro-Biopsy Cassette in QUICKLOAD™ Sleeves
 - M535-2SL – White
 - M535-3SL – Pink
 - M535-4SL – Green
 - M535-5SL – Yellow
 - M535-6SL – Blue
 - M535-7SL – Peach
 - M535-8SL – Tan
 - M535-9SL – Gray
 - M535-10SL – Lilac
 - M535-11SL – Orange
 - M535-12SL – Aqua
 - **M535AT** - SUPERFLO™ Tissue Cassette with Additive in QUICKLOAD™ Stack (Taped)
 - M535-2AT – White Laser Printer
 - M535-3AT – Pink Laser Printer
 - M535-4AT – Green Laser Printer
 - M535-5AT – Yellow Laser Printer
 - M535-6AT – Blue Laser Printer
 - M535-7AT – Peach Laser Printer
 - M535-8AT – Tan Laser Printer
 - M535-9AT – Gray Laser Printer
 - M535-10AT – Lilac Laser Printer
 - M535-11AT – Orange Laser Printer
 - M535-12AT – Aqua Laser Printer
- SWINGSETTE™ Series**
- **M515** - SWINGSETTE™ Tissue Cassettes
 - M515-2 – White
 - M515-3 – Pink
 - M515-4 – Green
 - M515-5 – Yellow
 - M515-6 – Blue
 - M515-7 – Peach
 - M515-8 – Tan
 - M515-9 – Gray
 - M515-10 – Lilac
 - M515-11 – Orange
 - M515-12 – Aqua
 - **M516** - SWINGSETTE™ Biopsy Cassettes
 - M516-2 – White
 - M516-3 – Pink
 - M516-4 – Green
 - M516-5 – Yellow
 - M516-6 – Blue
 - M516-7 – Peach
 - M516-8 – Tan
 - M516-9 – Gray
 - M516-10 – Lilac
 - M516-11 – Orange
 - M516-12 – Aqua
 - **M517** - SWINGSETTE™ Tissue Cassettes
 - M517-2 – White Laser Printer
 - M517-2A – White Laser Printer
 - M517-3 – Pink Laser Printer
 - M517-3A – Pink Laser Printer
 - M517-4 – Green Laser Printer
 - M517-4A – Green Laser Printer
 - M517-5 – Yellow Laser Printer
 - M517-5A – Yellow Laser Printer
 - M517-6 – Blue Laser Printer
 - M517-6A – Blue Laser Printer
 - M517-7 – Peach Laser Printer
 - M517-7A – Peach Laser Printer
 - M517-8 – Tan Laser Printer
 - M517-8A – Tan Laser Printer
 - M517-9 – Gray Laser Printer
 - M517-9A – Gray Laser Printer
 - M517-10 – Lilac Laser Printer
 - M517-10A – Lilac Laser Printer
 - M517-11 – Orange Laser Printer
 - M517-11A – Orange Laser Printer
 - M517-12 – Aqua Laser Printer
 - M517-12A – Aqua Laser Printer
 - **M517SL** - SWINGSETTE™ Tissue Cassettes in QUICKLOAD™ Sleeves
 - M517-2SL – White
 - M517-3SL – Pink
 - M517-4SL – Green
 - M517-5SL – Yellow
 - M517-6SL – Blue
 - M517-7SL – Peach
 - M517-8SL – Tan
 - M517-9SL – Gray
 - M517-10SL – Lilac
 - M517-11SL – Orange
 - M517-12SL – Aqua
 - **M517T** - SWINGSETTE™ Tissue Cassettes in QUICKLOAD™ Stack (Taped)
 - M517-2T – White Laser Printer
 - M517-2AT – White Laser Printer
 - M517-3T – Pink Laser Printer
 - M517-3AT – Pink Laser Printer
 - M517-4T – Green Laser Printer
 - M517-4AT – Green Laser Printer
 - M517-5T – Yellow Laser Printer
 - M517-5AT – Yellow Laser Printer
 - M517-6T – Blue Laser Printer
 - M517-6AT – Blue Laser Printer
 - M517-7T – Peach Laser Printer
 - M517-7AT – Peach Laser Printer
 - M517-8T – Tan Laser Printer
 - M517-8AT – Tan Laser Printer
 - M517-9T – Gray Laser Printer
 - M517-9AT – Gray Laser Printer
 - M517-10T – Lilac Laser Printer
 - M517-10AT – Lilac Laser Printer
 - M517-11T – Orange Laser Printer
 - M517-11AT – Orange Laser Printer
 - M517-12T – Aqua Laser Printer
 - M517-12AT – Aqua Laser Printer
 - **M518** - SWINGSETTE™ Biopsy Cassettes
 - M518-2 – White Laser Printer
 - M518-2A – White Laser Printer
 - M518-3 – Pink Laser Printer
 - M518-3A – Pink Laser Printer
 - M518-4 – Green Laser Printer
 - M518-4A – Green Laser Printer
 - M518-5 – Yellow Laser Printer
 - M518-5A – Yellow Laser Printer
 - M518-6 – Blue Laser Printer
 - M518-6A – Blue Laser Printer
 - M518-7 – Peach Laser Printer

- M518-7A – Peach Laser Printer
- M518-8 – Tan
- M518-8A – Tan Laser Printer
- M518-9 – Gray
- M518-9A – Gray Laser Printer
- M518-10 – Lilac
- M518-10A – Lilac Laser Printer
- M518-11 – Orange
- M518-11A – Orange Laser Printer
- M518-12 – Aqua
- M518-12A – Aqua Laser Printer
- **M518SL - SWINGSETTE™** Biopsy Cassettes in QUICKLOAD™ Sleeves
 - M518-2SL – White
 - M518-3SL – Pink
 - M518-4SL – Green
 - M518-5SL – Yellow
 - M518-6SL – Blue
 - M518-7SL – Peach
 - M518-8SL – Tan
 - M518-9SL – Gray
 - M518-10SL – Lilac
 - M518-11SL – Orange
 - M518-12SL – Aqua
- **M518T - SWINGSETTE™** Biopsy Cassettes in QUICKLOAD™ Stack (Taped)
 - M518-2T – White
 - M518-2AT – White Laser Printer
 - M518-3T – Pink
 - M518-3AT – Pink Laser Printer
 - M518-4T – Green
 - M518-4AT – Green Laser Printer
 - M518-5T – Yellow
 - M518-5AT – Yellow Laser Printer
 - M518-6T – Blue
 - M518-6AT – Blue Laser Printer
 - M518-7T – Peach
 - M518-7AT – Peach Laser Printer
 - M518-8T – Tan
 - M518-8AT – Tan Laser Printer
 - M518-9T – Gray
 - M518-9AT – Gray Laser Printer
- M518-10T – Lilac
- M518-10AT – Lilac Laser Printer
- M518-11T – Orange
- M518-11AT – Orange Laser Printer
- M518-12T – Aqua
- M518-12AT – Aqua Laser Printer
- **MACROSETTE® Series**
 - **M512 - MACROSETTE®** Processing Embedding Cassettes
 - M512 – White
- **MICROSETTE® Series**
 - **M502 - MICROSETTE®** Biopsy Cassettes
 - M502-2 – White
 - M502-2A – White Laser Printer
 - M502-3 – Pink
 - M502-3A – Pink Laser Printer
 - M502-4 – Green
 - M502-4A – Green Laser Printer
 - M502-5 – Yellow
 - M502-5A – Yellow Laser Printer
 - M502-6 – Blue
 - M502-6A – Blue Laser Printer
 - M502-7 – Peach
 - M502-7A – Peach Laser Printer
 - M502-8 – Tan
 - M502-8A – Tan Laser Printer
 - M502-9 – Gray
 - M502-9A – Gray Laser Printer
 - M502-10 – Lilac
 - M502-10A – Lilac Laser Printer
 - M502-11 – Orange
 - M502-11A – Orange Laser Printer
 - M502-12 – Aqua
 - M502-12A – Aqua Laser Printer
 - **M502AT - MICROSETTE®** Biopsy Cassettes with Additive in QUICKLOAD™ Stack (Taped)
 - M502-2AT – White Laser Printer
 - M502-3AT – Pink Laser Printer
- **M503 - MICROSETTE®** Biopsy Cassettes with 6 Compartments
 - M503-2 – White
 - M503-3 – Pink
 - M503-4 – Green
 - M503-5 – Yellow
 - M503-6 – Blue
 - M503-7 – Peach
 - M503-8 – Tan
 - M503-9 – Gray
 - M503-10 – Lilac
 - M503-11 – Orange
 - M503-12 – Aqua
- M502-4AT – Green Laser Printer
- M502-5AT – Yellow Laser Printer
- M502-6AT – Blue Laser Printer
- M502-7AT – Peach Laser Printer
- M502-8AT – Tan Laser Printer
- M502-9AT – Gray Laser Printer
- M502-10AT – Lilac Laser Printer
- M502-11AT – Orange Laser Printer
- M502-12AT – Aqua Laser Printer
- **M480 Series - Embedding Cassettes**
 - **M480 - Embedding Cassettes**
 - M480-2 – White
 - M480-2A – White Laser Printer
 - M480-3 – Pink
 - M480-3A – Pink Laser Printer
 - M480-4 – Green
 - M480-4A – Green Laser Printer
 - M480-5 – Yellow
 - M480-5A – Yellow Laser Printer
 - M480-6 – Blue
 - M480-6A – Blue Laser Printer
 - M480-7 – Peach
 - M480-7A – Peach Laser Printer
 - M480-8 – Tan
 - M480-8A – Tan Laser Printer
 - M480-9 – Gray
 - M480-9A – Gray Laser Printer

- M480-10 – Lilac
 - M480-10A – Lilac Laser Printer
 - M480-11 – Orange
 - M480-11A – Orange Laser Printer
 - M480-12 – Aqua
 - M480-12A – Aqua Laser Printer
 - **M480T** – Embedding Cassettes in QUICKLOAD™ Stack (Taped)
 - M480-2T – White
 - M480-2AT – White Laser Printer
 - M480-3T – Pink
 - M480-3AT – Pink Laser Printer
 - M480-4T – Green
 - M480-4AT – Green Laser Printer
 - M480-5T – Yellow
 - M480-5AT – Yellow Laser Printer
 - M480-6T – Blue
 - M480-6AT – Blue Laser Printer
 - M480-7T – Peach
 - M480-7AT – Peach Laser Printer
 - M480-8T – Tan
 - M480-8AT – Tan Laser Printer
 - M480-9T – Gray
 - M480-9AT – Gray Laser Printer
 - M480-10T – Lilac
 - M480-10AT – Lilac Laser Printer
 - M480-11T – Orange
 - M480-11AT – Orange Laser Printer
 - M480-12T – Aqua
 - M480-12AT – Aqua Laser Printer
 - **M480SL** - Embedding Cassettes in QUICKLOAD™ Sleeves
 - M480-2SL – White
 - M480-3SL – Pink
 - M480-4SL – Green
 - M480-5SL – Yellow
 - M480-6SL – Blue
 - M480-7SL – Peach
 - M480-8SL – Tan
 - M480-9SL – Gray
 - M480-10SL – Lilac
 - M480-11SL – Orange
 - M480-12SL – Aqua
 - **M481** - Metal Lid for M480 Cassettes
- Histology Accessories**
- **M460** - Embedding Rings
 - M460 – White
 - M460-3 – Pink
 - M460-4 – Green
 - M460-5 – Yellow
 - M460-6 – Blue
 - **M470** - Tissue Capsules
 - **M471** - Tissue Capsules
- **M474** - Stainless Steel Base Molds
 - M474-1 – 7 x 7 x 5 mm
 - M474-2 – 15 x 15 x 5 mm
 - M474-3 – 24 x 24 x 5 mm
 - M474-4 – 30 x 24 x 5 mm
 - M474-5 – 37 x 24 x 5 mm
 - **M475** - Disposable Base Molds
 - M475-1 – 7 x 7 x 5 mm
 - M475-2 – 15 x 15 x 5 mm
 - M475-3 – 24 x 24 x 5 mm
 - M475-4 – 30 x 24 x 5 mm
 - M475-5 – 37 x 24 x 5 mm
 - **M475-10** - Deep Base Mold, Disposable
 - **M476** - Biopsy Foam Pads
 - M476-1 – Cassettes
 - M476-2 – Small Capsules
 - M476-3 – Large Capsules
 - M476-4 – Slimsette
 - M476-5 – Histosette
 - **M478** - Biopsy Bags
 - M478-1 – 30 x 50 mm
 - M478-2 – 45 x 75 mm
 - M478-3 – 75 x 95 mm

Appendix II: Revision History Table

Rev	Revision Date	Description of Change	Originator
1.0	03/05/2022	Initial Release – Declaration of Conformity	Annette Roy
1.1	19/09/2024	Addition of Appendix I: Addendum to the Original Declaration of Conformity	Rossano Morgan-Cavallaro
1.2	01/10/2024	The Addition of Appendix I – List of Trade Names and Product Reference Codes and it's alignment with the Product Listing Addendum.	Rossano Morgan-Cavallaro
1.3	16/04/2025	Reformatted to EMERGO template with Revision History Table. No change to Declaration content.	Rossano Morgan-Cavallaro
1.4	20/08/2025	Followed Emergo's expert recommendation to remove "Used by laboratory professionals only" from the Intended Purpose Statement. Also added Cassettes with Additives to applicable product series in Appendix I, now including variants optimized for laser printing. No change to form or function.	Rossano Morgan-Cavallaro