

TECHNICAL DATA SHEET
PRODUCT DESCRIPTION

PRODUCT CODE	147C
GENERAL DESCRIPTION	VIRO-SWAB T/LIQ.PLAS.
DETAILED DESCRIPTION	10 boxes containing 25 film-film pouches each containing 1 polyester tip applicator with fuchsia cap and a tube with white cap and sponge with 1.2 ml of medium for viral infections (Chlamydia)

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
VIRAL MEDIUM	Sodium chloride Potassium chloride Monopotassium phosphate Sodium thioglycollate Calcium Chloride Magnesium chloride Disodium phosphate Distilled Water

APPLICATOR	PLASTIC APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH: 144.0 mm \pm 1.0 mm PLASTIC TIP DIAMETER: 2.50 mm \pm 0.20 mm FIBRE TIP LENGTH: 15.0 mm \pm 2.0 mm FIBRE TIP DIAMETER: 4.50 mm \pm 0.50 mm TOTAL LENGTH: 147.0 mm \pm 1.0 mm
TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 154.50 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.2 mm EXTERNAL DIAMETER: 12.0 mm \pm 0.5 mm



TECHNICAL DATA SHEET

SPONGE	SPONGE	DIMENSIONAL SPECIFICATIONS
		LENGTH: 37.50 mm + 2.0 mm -1.50 mm WIDTH: 10.0 mm ± 0.5 mm THICKNESS: 5.0 mm ± 0.5 mm
CAP	WHITE CAP W/O RECEPTACLE	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.0 mm ± 0.10 mm CAP LENGTH: 32.80 mm ± 0.20 mm INSERTION DIAMETER: 9.00 mm +0.05 mm
CAP	CAP PINK BIG HOLE	DIMENSIONAL SPECIFICATIONS
		EXTERNAL CAP DIAMETER: 12.00 mm ± 0.20 mm CAP LENGTH: 32.90 mm ± 0.50 mm INSERTION DIAMETER: 9.30 mm + 0.20 mm

PACKAGING SPECIFICATION

INNER BOX QTY	25
OUTER BOX QTY	250
PACKAGING	FILM-FILM POUCH

**TECHNICAL DATA SHEET****IMPORTANT INFO**

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	15 MONTHS
MEDIUM pH	7.2 ± 0.3