

TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	148C
GENERAL DESCRIPTION	VIRO-SWAB T/LIQ.ALL.
DETAILED DESCRIPTION	10 boxes containing 25 film-film pouches each containing an aluminium applicator with polyester tip with fuchsia cap and a tube with white cap without receptacle and sponge with 1.2 ml of medium for viral infections (Chlamydia)

PRODUCT COMPONENTS

MEDIUM		COMPOSITION
VIRAL MEDIUM		Sodium chloride Potassium chloride Monopotassium phosphate Sodium thioglycollate Calcium Chloride Magnesium chloride Disodium phosphate Distilled Water
APPLICATOR	ALUMINIUM APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH: 150.0 mm \pm 1.0 mm TOTAL LENGTH: 152.0 mm \pm 1.0 mm STICK DIAMETER: 1.00 mm \pm 0.10 mm TIP DIAMETER: 3.00 mm \pm 0.50 mm TIP LENGTH: 15.00 mm \pm 2.00 mm
CAP	FUCHSIA CAP	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.20 mm \pm 0.10 mm CAP LENGTH: 32.50 mm \pm 0.20 mm INSERTION DIAMETER: 9.30 mm \pm 0.05 mm



TECHNICAL DATA SHEET

TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 154.50 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.2 mm EXTERNAL DIAMETER: 12.0 mm \pm 0.5 mm
SPONGE	SPONGE	DIMENSIONAL SPECIFICATIONS
		LENGTH: 37.50 mm + 2.0 mm -1.50 mm WIDTH: 10.0 mm \pm 0.5 mm THICKNESS: 5.0 mm \pm 0.5 mm
CAP	WHITE CAP W/O RECEPTACLE	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.0 mm \pm 0.10 mm CAP LENGTH: 32.80 mm \pm 0.20 mm INSERTION DIAMETER: 9.00 mm +0.05 mm

PACKAGING SPECIFICATION

INNER BOX QTY	25
OUTER BOX QTY	250
PACKAGING	FILM-FILM POUCH

**TECHNICAL DATA SHEET****IMPORTANT INFO**

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	15 MONTHS
MEDIUM pH	7.2 ± 0.3