

TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	149C
GENERAL DESCRIPTION	VIRO-SWAB T/LIQ. TWISTED APP.
DETAILED DESCRIPTION	10 boxes containing 25 film-film pouches containing 1 tube with sponge, 1.2mL of medium for VIRAL / Chlamydia and white cap without receptacle and 1 tube with aluminum applicator with polyester tip and fuchsia cap.

PRODUCT COMPONENTS

MEDIUM		COMPOSITION
VIRAL MEDIUM		Sodium chloride Potassium chloride Monopotassium phosphate Sodium thioglycollate Calcium Chloride Magnesium chloride Disodium phosphate Distilled Water
APPLICATOR	ALUMINIUM APPLICATOR	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 150.00 mm \pm 1 mm STICK DIAMETER HANDLE PART: 1.6 mm \pm 0.10 mm TIP DIAMETER: 2.80 mm \pm 0.50 mm TIP LENGTH: 16.00 mm \pm 5 mm
CAP	FUCHSIA CAP	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.20 mm \pm 0.10 mm CAP LENGTH: 32.50 mm \pm 0.20 mm INSERTION DIAMETER: 9.30 mm \pm 0.05 mm



TECHNICAL DATA SHEET

TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 155.00 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.5 mm EXTERNAL TUBE DIAMETER: 12.0 mm \pm 0.5 mm
TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 154.50 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.2 mm EXETERNAL DIAMETER: 12.0 mm \pm 0.5 mm
SPONGE	SPONGE	DIMENSIONAL SPECIFICATIONS
		LENGTH: 37.50 mm + 2.0 mm -1.50 mm WIDTH: 10.0 mm \pm 0.5 mm THICKNESS: 5.0 mm \pm 0.5 mm
CAP	WHITE CAP W/O RECEPTACLE	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 12.0 mm \pm 0.10 mm CAP LENGTH: 32.80 mm \pm 0.20 mm INSERTION DIAMETER: 9.00 mm +0.05 mm

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

INNER BOX QTY	25
OUTER BOX QTY	250
PACKAGING	FILM-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	15 MONTHS
MEDIUM pH	7.2 ± 0.3