



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

|                             |                                                                                                                                                                                                                        |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>PRODUCT CODE</b>         | 149C                                                                                                                                                                                                                   |
| <b>GENERAL DESCRIPTION</b>  | VIRO-SWAB T/LIQ. TWISTED APP.                                                                                                                                                                                          |
| <b>DETAILED DESCRIPTION</b> | 10 boxes containing 25 film-film pouches containing 1 tube with sponge, 1.2mL of medium for VIRAL / Chlamydia and white cap without receptacle and 1 tube with aluminum applicator with polyester tip and fuchsia cap. |

### **PRODUCT COMPONENTS**

| <b>MEDIUM</b>     |                             | <b>COMPOSITION</b>                                                                                                                                                           |
|-------------------|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| VIRAL MEDIUM      |                             | Sodium chloride<br>Potassium chloride<br>Monopotassium phosphate<br>Sodium thioglycollate<br>Calcium Chloride<br>Magnesium chloride<br>Disodium phosphate<br>Distilled Water |
| <b>APPLICATOR</b> | <b>ALUMINIUM APPLICATOR</b> | <b>DIMENSIONAL SPECIFICATIONS</b>                                                                                                                                            |
|                   |                             | TOTAL LENGTH: 150.00 mm $\pm$ 1 mm<br>STICK DIAMETER HANDLE PART: 1.6 mm $\pm$ 0.10 mm<br>TIP DIAMETER: 2.80 mm $\pm$ 0.50 mm<br>TIP LENGTH: 16.00 mm $\pm$ 5 mm             |
| <b>CAP</b>        | <b>FUCHSIA CAP</b>          | <b>DIMENSIONAL SPECIFICATIONS</b>                                                                                                                                            |
|                   |                             | EXTERNAL DIAMETER: 12.20 mm $\pm$ 0.10 mm<br>CAP LENGTH: 32.50 mm $\pm$ 0.20 mm<br>INSERTION DIAMETER: 9.30 mm $\pm$ 0.05 mm                                                 |



## TECHNICAL DATA SHEET

|                       |                             |                                                                                                                                 |
|-----------------------|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>TRANSPORT TUBE</b> | TUBE                        | <b>DIMENSIONAL SPECIFICATIONS</b>                                                                                               |
|                       |                             | TUBE LENGTH: 155.00 mm $\pm$ 1.00 mm<br>INSERTION DIAMETER: 9.0 mm $\pm$ 0.5 mm<br>EXTERNAL TUBE DIAMETER: 12.0 mm $\pm$ 0.5 mm |
| <b>TRANSPORT TUBE</b> | TUBE                        | <b>DIMENSIONAL SPECIFICATIONS</b>                                                                                               |
|                       |                             | TUBE LENGTH: 154.50 mm $\pm$ 1.00 mm<br>INSERTION DIAMETER: 9.0 mm $\pm$ 0.2 mm<br>EXTERNAL DIAMETER: 12.0 mm $\pm$ 0.5 mm      |
| <b>SPONGE</b>         | SPONGE                      | <b>DIMENSIONAL SPECIFICATIONS</b>                                                                                               |
|                       |                             | LENGTH: 37.50 mm + 2.0 mm -1.50 mm<br>WIDTH: 10.0 mm $\pm$ 0.5 mm<br>THICKNESS: 5.0 mm $\pm$ 0.5 mm                             |
| <b>CAP</b>            | WHITE CAP W/O<br>RECEPTACLE | <b>DIMENSIONAL SPECIFICATIONS</b>                                                                                               |
|                       |                             | EXTERNAL DIAMETER: 12.0 mm $\pm$ 0.10 mm<br>CAP LENGTH: 32.80 mm $\pm$ 0.20 mm<br>INSERTION DIAMETER: 9.00 mm +0.05 mm          |

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

|                      |                 |
|----------------------|-----------------|
| <b>INNER BOX QTY</b> | 25              |
| <b>OUTER BOX QTY</b> | 250             |
| <b>PACKAGING</b>     | FILM-FILM POUCH |

**IMPORTANT INFO**

|                             |                      |
|-----------------------------|----------------------|
| <b>STERILIZATION METHOD</b> | GAMMA IRRADIATION    |
| <b>STORAGE TEMPERATURE</b>  | 5/25°C-DO NOT FREEZE |
| <b>SHELF LIFE</b>           | 15 MONTHS            |
| <b>MEDIUM pH</b>            | 7.2 ± 0.3            |