



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	8177CS10H
GENERAL DESCRIPTION	LOOPS 10UL HARD 2000PKG
DETAILED DESCRIPTION	1 outer box contain 4 inner boxes each composed of 50 paper-film pouches each containing 10 10uL hard loops.

PRODUCT COMPONENTS

		SPECIFICATIONS
LOOP	LOOP 10 UL	TOTAL LENGTH: 201.00 mm ± 0.50 mm INSIDE DIAMETER RING: 3.80 mm ± 0.10 mm OUTER DIAMETER RING: 6.50 mm ± 0.10 mm THICK RING: 1.20 mm ± 0.10 mm MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	500
OUTER BOX QTY	2000
PACKAGING	PAPER-FILM POUCH
ENVELOPE QTY	10

IMPORTANT INFO

STERILIZATION METHOD	IONIZING IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.02-05/05/2022) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.