

Description


- Whirl-Pak® Write On Sterilized Bags

Part number

- 052920

Pieces by sales :

500 Bags



Description

Whirl-Pak® Write-On Sterilized Sampling Bags provide a secure, contaminant free flexible container that ensures dependable analysis results. You can count on Whirl-Pak® products as a proven, economical and efficient way to collect, store and transport samples.

Technical data

- Wire Style: round
- Volume – fl oz: 55
- Volume – mL: 1,627
- Dimensions – inch (WxL): 7,5 x 12
- Dimensions – cm (WxL): 19 x 30
- Film Thickness: 4

