

Brand: Honeywell Fluka
 Product: 06203

Sodium hydroxide pellets, puriss., meets analytical specification of Ph. Eur., BP, NF, E 524

Parameter	Specification	Units
assay	98.0-100.5	%
assay of Na ₂ CO ₃	max. 0.5	%
aluminium (Al)	max. 0.001	%
arsenic (As)	max. 2	ppm
calcium (Ca)	max. 0.001	%
iron (Fe)	max. 10	ppm
mercury (Hg)	max. 1	ppm
potassium (K)	max. 0.1	%
lead (Pb)	max. 0.5	ppm
heavy metals (as Pb)	max. 20	ppm
chloride (Cl)	max. 100	ppm
total N	max. 0.0005	%
phosphate (PO ₄)	max. 0.001	%
silicate (as SiO ₂)	max. 0.01	%
sulfate (SO ₄)	max. 100	ppm
water insoluble matter + org. comp.	complying	
residual solvents	complying	
appearance of the solution	complying	