

# hybex™

## Media Storage Bottles



**Borosilicate  
Glass**

hybex™ reusable media bottles are manufactured from high quality, borosilicate 3.3 glass. With their superior chemical resistance, these bottles are ideal for storage of reagents, culture media, biological fluids and a variety of other aqueous and nonaqueous solutions.

A low coefficient of thermal expansion (3.3), provides temperature resistance that is far superior to that of standard lab glass (or soda lime glass). Hybex glass can withstand temperatures from -70 to 500°C, permitting repeated use in freezers, microwaves, autoclaves and dry heat sterilization chambers.

All hybex bottles have permanent, white enamel graduations for volume measurement and a white identification box for marking or coding.

All hybex media bottles include a drip-free sealing ring as well as a standard blue polypropylene (GL45) sealing cap. Bottles are also available with green, purple, white and yellow caps, or in an assortment pack containing two of each color.

— **Borosilicate 3.3 hybex glass**

— **Autoclavable caps (GL45) & drip rings**

— **Permanent enamel graduation markings**

— **Large area for bottle marking/identification**

— **Autoclavable, reuseable & breakage resistant**

— **Five cap colors**

### Technical Data

|                  |                                                  |
|------------------|--------------------------------------------------|
| Material:        | hybex™ borosilicate glass (3.3)                  |
| Cap Thread Size: | GL45*                                            |
| Temp. Range:     | -70°C to 500°C                                   |
| Sterilization:   | Autoclave or dry heat                            |
| Cap Material:    | Polypropylene (140°C)                            |
| Cap Colors:      | Blue, green, purple<br>yellow, white or assorted |

\*50ml bottles are only supplied with blue GL32 sealing caps

### Ordering Information

|             |                                                                                                                                                     |
|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| B3000-50    | Media storage bottle, 50ml with standard (GL32) blue cap, 10/pk.                                                                                    |
| B3000-100*  | Media storage bottle, 100ml with standard (GL45) blue cap, 10/pk.                                                                                   |
| B3000-250*  | Media storage bottle, 250ml with standard (GL45) blue cap, 10/pk.                                                                                   |
| B3000-500*  | Media storage bottle, 500ml with standard (GL45) blue cap, 10/pk.                                                                                   |
| B3000-1000* | Media storage bottle, 1L with standard (GL45) blue cap, 10/pk.                                                                                      |
| B3000-2000* | Media storage bottle, 2L with standard (GL45) blue cap, 5/pk.                                                                                       |
| B3000-SP*   | Media storage bottle, <b>Starter Pack</b> with standard (GL45) blue cap (Starter Pack of ten bottles includes 2x100ml, 3x250ml, 3x500ml & 2x1000ml) |

|            |                                                      |
|------------|------------------------------------------------------|
| B3000-CAP* | High temperature red caps, up to 180°C (GL45), 10/pk |
| B3000-HTC  | Replacement cap, standard blue (GL45), 10/pk.        |
| B3000-RIN  | Replacement sealing ring, (GL45), 10/pk.             |
| B3000-CAP2 | Replacement cap, standard blue (GL32), 10/pk.        |
| B3000-RIN2 | Replacement sealing ring, (GL32), 10/pk.             |

\*Please add **(-G)** for green caps, **(-P)** for purple caps, **(-Y)** for yellow caps, **(-W)** for white caps or **(-AS)** for assorted caps (2 of ea. color)

**Benchmark**  
Scientific

**Benchmark Scientific Inc.**  
PO Box 709

Edison, NJ 08818

Phone: 908-769-5555

Fax: 908-222-1864

Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)

Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)