



Drummond Scientific Company

500 PARKWAY, BROOMALL, PENNSYLVANIA 19008-4203

PHONE : 610 – 353 – 0200

FAX 610 – 353 – 6204

www.drummondsci.com

Drummond Scientific 5ul Calibrated Pipets Specification Sheet

Catalog # 2-000-005

Color code: White

5ul Line Location: 2.165”

Length: 5.000”

OD: .0430”

ID: .01340”

Glass type: Borosilicate

Packaging: 250 pieces per container

Aspirator Tube included

DOMINIQUE DUTSCHER SAS