



Drummond Scientific Company Technical Data Sheet

**Product # 1-000-0010**

**1.0 ul Microcaps**

**Dimensions**

Volume: 1.0 ul (+/- 1.00%)

Length: 1.260" +/- .003" (32mm +/- .076mm)

Glass O.D. – 0.260" +/- .0020" (.660mm +/- .051mm)

Glass I.D. - .00785" +/- .00004" (.1994mm +/- .0010mm)

Glass type: Borosilicate

**Packaging**

Container: Printed Drummond Flex Tube

Vial Size: 9/16" dia. x 2.00" L

Foam Plug: 1/2" dia. x 3/4" L

Cap Size: 15mm w/ .040 hole

Quantity per vial: 100 glass pipettes.

Packaged with Bulb Assembly

Material listed herein is produced in accordance with ASTM Standard Specification for Disposable Glass Micropipets number E 672-87. Standards used in manufacture are traceable to the National Institute of Standards and Technology.

Alexis Piccirilli

Director of Quality

Drummond Scientific Company