

**Petri dish, 55mm diameter, without rim, s/15, aseptic, Clearline x1200**

- Description: Polystyrene Petri dishes for microbiology, bacteriology
- Made of crystal polystyrene, superior clarity
- Manufactured under aseptic conditions, ISO Class 7
- Available aseptic or sterile (E-beam irradiated, 10 kGy)
- Sub-packaged in bags
- Suitable for food contact, according to EC Regulation No. 1935/2004
- Pouring temperature 60°C
  
- Model Ø 55
- Lugs: None
- Diameter (mm) Internal/External : 55/58
- Height (mm) Internal/External : 13.1/14.2
- Weight (g) : 6.1
- Units per bag 15
- Sterility : Aseptic

**Units : 1200**