

BEAD RUPTOR 24 7 ML TUBE CARRIAGE KIT

Sku: 19-345-007



GENERAL INFO

Bead Ruptor 24 7 mL Tube Carriage Kit

7 mL tube carriage holds 12 x 7 mL screw cap tubes

Ordering Information: Includes one 12 x 7 mL screw cap tube carriage and one lock ring.

SPECIFICATIONS

Capacity:	12 x 7 mL skirted screw cap tubes
Volume Range:	0.5 mL - 5.0 mL
Compatibility:	Requires Omni 7 mL tubes
Force Limit:	Maximum force of 6 m/s